M-PM[®] Disc

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- | | |
|--------------------------------------|------------------------------------|
| (en) Instructions for use | (lt) Naudojimo instrukcija |
| (fr) Notice d'utilisation | (lv) Lietošanas informācija |
| (it) Istruzioni per l'uso | (et) Kasutusjuhend |
| (es) Instrucciones de uso | (hu) Használati utasítás |
| (pt) Instruções de utilização | (cs) Návod k použití |
| (pl) Instrukcja użycia | (sl) Navodila za uporabo |
| (nl) Gebruiksaanwijzing | (sk) Návod na používanie |
| (sv) Bruksanvisning | (bg) Инструкции за употреба |
| (da) Brugsanvisningen | (el) Οδηγίες χρήσης |
| (no) Bruksanvisningen | (tr) Kullanım kılavuzu |
| (fi) Käyttöohje | |

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M-PM® Disc, PMMA milling blank without fillers and fibres

tooth colour, pink colour, clear

(en) Instructions for use, please read carefully!

Intended purpose

- **M-PM-Disc, tooth colour**, for functional restoration by replacement of lost tooth substance in permanent dentition in the form of temporary crowns and bridges in the anterior and posterior regions.

Indication

- Loss of hard tooth substance and resulting aesthetic impairment

- **M-PM® Disc, pink colour**, replacement of lost soft tissue in the oral region with a denture base

Indication

- Total or partial tooth loss

- **M-PM® Disc, clear**, removal of overburden or incorrect burden of hard tooth substance, the periodontal apparatus and the mandibular joints with splints (occlusal splints) / fabrication of patient-customised drilling templates for the correct positioning of implants

Indication

- Physiological dysfunction leading to abrasion of hard tooth substance

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Intended users

Dental professionals

Contraindication M-PM® Disc, tooth colour, pink colour and clear

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Composition

M-PM® Disc, tooth colour, pink colour

- Polymethylmethacrylate (PMMA) and cross-linked polymers based on methacrylic acid esters
- Colourants
- Residual peroxide as dibenzoyl peroxide
- Methylmethacrylate (MMA) may be contained as residual monomer up to max. 1%.

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General information on M-PM® Disc constructions (tooth coloured, pink, clear)

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Special information on M-PM® Disc constructions, tooth colour

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Veneering of PMMA frameworks, crowns, bridges (M-PM® Disc, tooth colour)

- Milled PMMA frameworks, crowns and bridges may be coated with conventional veneering materials. Materials based on MMA are preferable here.
- When veneering with composites the connection with the PMMA must always be made according to the respective manufacturer's specifications.

Trimming

- Cross-cut carbide burs suitable for resin or suitable cutting discs should be used to rework and separate frameworks, crowns, bridges, denture bases, bite splints and drilling templates from the blank.
- To prevent inaccurate fitting after fabricating the construction, avoid excessive heat development during trimming and polishing.

Restoration M-PM® Disc, tooth colour, pink colour, clear

Frameworks, crowns, bridges, denture bases, bite splints and drilling templates, which have been fabricated from the M-PM® Disc can be restored and supplemented at any time with all standard MMA-based cold-curing polymers, such as Combipress N/LM, Weropress® or artConcept®. Prior conditioning of the surface is necessary, preferably with artConnect or a similar MMA-based bonding agent.

Fixation in the mouth (M-PM® Disc, tooth colour)

- | | |
|-------------------------------|---|
| for the temporary restoration | for long-term temporary restorations |
| - temporary cements | - glass ionomer cements without bonding |

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Warnings

Processing PMMA frameworks may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system on your milling machine and at the workplace is working properly for individual post-processing and that you are using personal protective equipment.

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en

can

Instructions for use, please read carefully!

Purpose of use

PMMA-Disc is a milling blank to fabricate

- acrylic crowns, bridges and framework for provisional restorations within the prosthetic treatment of anterior and posterior teeth.
- In constructing bridges the maximum number of units is 6, number of pontics 2. Cantilevered pontics: none.

Composition

- Polymethyl methacrylate (PMMA) and cross-linking polymers on basis of methacrylic acid esters
- Dyes
- Residual peroxide
- Methyl methacrylate (MMA) can be contained as residual monomer by up to 1 % max.

Instruction for the user

Construction

- In constructing bridges the maximum number of units is 6, number of pontics 2. Cantilevered pontics: none.
- The minimum wall thickness of 1 mm must not fall below.
- A span width of more than one span unit is contraindicated in bridge constructions. The connection area between the crown and bridge span must not fall below 9 mm².

- The user must decide in individual cases, according to the relevant indication, whether the PMMA - OMP-N® material is providing sufficient masticatory stability required for the restoration, when selecting a bonding system.

Veneering PMMA frameworks

- A milled PMMA framework can be coated with conventional veneering materials. It is recommended to use materials on PMMA basis.
- When using a composite veneering material, the bonding to the PMMA must always be performed according to the instructions of the respective manufacturer.

Finishing

- For refinishing and removing the framework from the blank, it is recommended to use helical-tooth cutters (hard metal) or suitable cutting discs.
- In order to avoid problems with the marginal fit after constructing the framework, it is recommended to avoid excessive heat during the preparation and polishing.

Repair

- Frameworks, whole crowns and bridges manufactured from the PMMA-Discs can be repaired and added on to at any time with any customary cold-cure polymer on MMA basis.

Options for oral attachment

For the provisional restoration

- temporary cements

For long-term temporaries

- glass ionomer cements without bonding

Storage

- No special measures required

Special instructions

- Indicate the LOT No. for every procedure requiring an identification of the material.
- All serious incidents occurring with users and / or patients in connection with this product must be reported to Merz Dental, as well as to the German Federal Institute for Drugs and Medical Devices (BfArM), or the respective competent authority of the member state in which the user and / or patient is resident.

Safety precautions

- Dust is created during processing, which may lead to irritations of the eyes, skin and airways. Therefore, always pay attention that the exhaust system of your milling machine is working perfectly.
- Further safety information can be obtained from the material safety data sheet.

Disposal

- The PMMA-Disc is water-insoluble, inactive, not harmful to the groundwater and can therefore be disposed via the normal household waste collection.

The product characteristics are based on compliance and observance of this instruction for use.

Date of information 2021-11

DENTAL

Merz Dental GmbH

Kieferweg 1 24321 Lütjenburg, Germany

Tel + 49 (0) 4381 / 403-0

Fax + 49 (0) 4381 / 403-403

www.merz-dental.de

EN ISO 13485

M-PM[®] Disc, PMMA milling blank without fillers and fibres

tooth color, pink, clear

usa Instructions for use, please read carefully!

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed dental practitioner and dental technician.

Intended use

- **M-PM[®]-Disc, tooth color**, for use as a milling blank in the fabrication of provisional crowns and bridges.
- For use as a milling blank in the fabrication of permanent removable dental prostheses.
- **M-PM[®] Disc, pink** for use as a milling blank for the fabrication of denture bases for removable dentures.
- **M-PM[®] Disc, clear**, indicated for the fabrication of removable bite splints such as mouthguards and nightguards.

Classification according to MDD 93/42/EEC and MDR (EU) 2017/745

Class IIa medical device for long-term contact with intact mucous membranes or temporary contact with injured mucous membranes.

Patient target groups

All genders and age groups

Intended users

Dental professionals

Contraindication M-PM[®] Disc, tooth color, pink and clear

- Hypersensitivity to acrylate compounds

Composition

M-PM® Disc, tooth color, pink

- Polymethylmethacrylate (PMMA) and cross-linked polymers based on methacrylic acid esters
- Colorants
- Residual peroxide as dibenzoyl peroxide
- Methylmethacrylate (MMA) may be contained as residual monomer up to max. 1%.

M-PM® Disc, clear

- Polymethylmethacrylate (PMMA) and cross-linked polymers based on methacrylic acid esters
- Residual peroxide as dibenzoyl peroxide
- Methylmethacrylate (MMA) may be contained as residual monomer up to max. 1%.

General information on M-PM® Disc constructions (tooth colored, pink, clear)

- Minimum wall thickness after polishing ≥ 1 mm

Special information on M-PM® Disc constructions, tooth color

- In the construction of bridges with one pontic, the connection cross-section between the crown and pontic must not be less than 9 mm².
- In the construction of bridges spanning up to two pontics, the connection cross-section between crowns and pontics must not be less than 12 mm².
- By selecting the respective fixation system according to the indication, the user has to decide case-by-case whether the OMP-N® (PMMA) material ensures sufficient chewing stability for the restoration.

Veneering of PMMA substructures (M-PM® Disc, tooth color)

- A milled PMMA substructure can be veneered with conventional appropriate materials. It is recommended to use material on a MMA basis.
- For composite veneers, the bond to the PMMA must always be done according to the instructions of the manufacturer.

Notes for the user

Trimming

- Cross-cut carbide burs suitable for resin or suitable cutting discs should be used to work and separate crowns, bridges, denture bases, bite splints and drilling templates from the blank.
- To prevent inaccurate fitting after fabricating the construction, avoid excessive heat development during trimming and polishing.

Restoration M-PM® Disc, tooth color, pink color, clear

Constructions produced from the M-PM® Disc can be restored, supplemented or lined at any time with MMA-based cold-curing polymers, such as Weropress® in the matching shade. Following roughening of the surface (sanding / blasting / grinding), further conditioning is recommended with an MMA-based bonding agent.

Fixation in the mouth (M-PM® Disc, tooth color)

for the temporary restoration

for long-term temporary restorations

- temporary cements

- glass ionomer cements without bonding

Wear period for fixed dentures

- up to 3 years

Storage

- No special measures required. See label for expiration date.

Notes

- Changes to the construction through use, e.g. due to fissures, cracks, clouding or significant signs of abrasion, necessitate a repair or refabrication. Monitoring should be carried out with regular assessment of the denture at the follow-up presentation in the dental surgery.
- Record the LOT number with every process that requires identification of the material.
- To ensure traceability and assignment of the product data as well as the LOT numbers for the patient at any time, it is recommended to store the milling blank in the original packaging during two milling operations.
- All serious incidents occurring with users and / or patients in connection with this product must be reported to Merz Dental, as well as to the German Federal Institute for Drugs and Medical Devices (BfArM), or the respective competent authority of the member state in which the user and / or patient is resident.

Safety advice

Processing PMMA discs may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system on your milling machine and at the workplace is working properly for individual post-processing and that you are using personal protective equipment.

Warning: The fabrication of teeth for longterm removable dentures from M-PM® Disc (tooth color) must not be used in the preparation of bridges. For the preparation of bridges (maximum of two pontics) always refer to the IFU section on preparation of bridges from M-PM® Disc (tooth color).

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de.

Disposal

- Residues from processing the M-PM® Disc have to be disposed of as residual waste.

The product properties are based on compliance with and observation of these instructions for use.

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