artOpaque *PLUS* GUM

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- | | |
|------------------------------------|---------------------------------------|
| en Instructions for use | lt Naudojimo instrukcija |
| fr Notice d'utilisation | lv Lietošanas informācija |
| it Istruzioni per l'uso | et Kasutusjuhend |
| es Instrucciones de uso | hu Használati utasítás |
| pt Instruções de utilização | cs Návod k použití |
| pl Instrukcja użycia | sl Navodila za uporabo |
| nl Gebruiksaanwijzing | sk Návod na používanie |
| sv Bruksanvisning | bg Инструкции за употреба |
| da Brugsanvisningen | el Οδηγίες χρήσης |
| no Bruksanvisningen | tr Kullanım kılavuzu |
| fi Käyttöohje | uk Інструкція для застосування |

artOpaque^{PLUS} GUM

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artOpaque^{PLUS} GUM

Gingiva-coloured paste opaquer with integrated bonding component.

en Instructions for use, please read carefully!

Indication

- Total or partial tooth loss
- Loss of hard tooth substance

Intended use

Replacement of lost hard tooth substance by:

- Tooth-shaded fixture and for backing artVeneer®
- Veneering telescopic, attachment, bar and implant constructions with artVeneer®
- Tooth-shaded backing of trimmed acrylic teeth
- Fixing acrylic teeth to the model cast with a back protection plate or pontic
- Individualising acrylic teeth, artVeneer® and CAD/CAM restorations produced from M-PM® Disc and the artBloc® Temp
- Repairs

Contraindication

- Hypersensitivity to one of the components

Components of the artConcept® PLUS system

- artPrime PLUS - Metal bonding agent
- artPreOpaque PLUS - Priming opaquer
- artOpaque PLUS - Opaquer for covering framework
- artDentine PLUS/artEnamel PLUS Polymer - Tooth-shade PMMA acrylic
- artDentine PLUS Liquid NT (normal time) / artDentine PLUS Liquid LT (long time) - Monomer component for artDentine PLUS / artEnamel PLUS
- artConnect PLUS - Acrylic bonder
- artVeneer® - Tooth-shade PMMA anterior and posterior veneer

Composition

methacrylic acid ester liquid containing an activator, d,l camphor quinone

Processing

Processing of the model casts is as usual.

- Blast the veneer surface with aluminium oxide (50 µm, 2 bar) and blow off the activated surface with oil-free compressed air. The framework surfaces must be free of grease and abrasive residues.
- Apply artOpaque PLUS GUM thinly on the surface, as in washing, using a stiff shorthair brush and light cure. A non-covering opaquer layer has a negative impact on the shade effect. At least two coats are necessary in order to achieve complete coverage. Opaquer layers applied too thickly weaken the bond between framework and plastic base. Apply especially thinly in the region of the retentions.
- The fluidity of the artOpaque PLUS GUM can be optimised by stirring on a mixing block or similar device.
- The two available colour shades dark and light can be mixed freely for colour coordination.

- Polymerise each layer separately, extending the polymerisation time does not result in greater hardening.
- The processing time of the compounds is 1-3 min depending on the light conditions.

Curing unit artOpaque PLUS GUM

DC-XS / Uni-XS	Kulzer	90 sec
HeraFlash	Kulzer	90 sec
Labolight LV-II / III	GC	1 min
MPA 200	Dentsply	1 x HD
Solidilite EX / V	Shofu	1 min
Spektra 2000	Schütz Dental	3 min
Spektra 2000 fast	Schütz Dental	3 min
Spektra LED	Schütz Dental	60 sec
Spektramat	ivoclar vivadamat	2 min
Triad II	Dentsply	5 min

Storage

Store in a dry and cool area at 10 - 25 °C / 50 - 77 °F without direct exposure to light.

Disposal

The fully polymerized components of artConcept® PLUS are disposed of as residual waste. Disposal of non-polymerized product residues according to legal regulations. Do not allow entry into the sewage system or waterways.

Notes

- Close all containers tightly after use.
- Mixing with other opaquer fluids/powders is contraindicated.
- Do not use the material after the specified expiration date.
- To reduce the residual monomer content, the finished (sanded) but not polished polymerised workpiece should be stored in water at room temperature for about 48 hours before delivery.
- The residual monomer content after polymerisation of artDentine PLUS/artEnamel PLUS is < 1.8%.

Danger and safety instructions

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: Wash with plenty of water.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 Take off contaminated clothing and wash it before reuse.

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de.

The product properties are based on compliance with and observation of these instructions for use.

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Merz Dental GmbH