Baltic Denture System

Digital solution for dentures

BDKEY® Set

- | | |
|------------------------------------|----------------------------------|
| de Gebrauchsinformation | fi Käyttöohje |
| en Instructions for use | lt Informacinis lapelis |
| fr Mode d'emploi | lv Lietošanas informācija |
| it Istruzioni per l'uso | et Kasutusteave |
| es Instrucciones de uso | hu Használati tájékoztató |
| pt Instruções de utilização | cs Informace k použití |
| pl Instrukcja użycia | sl Navodila za uporabo |
| nl Gebruiksaanwijzing | sk Návod na použitie |
| sv Bruksanvisning | bg Информация за употреба |
| da Brugsanvisning | el Οδηγίες χρήσης |
| no Bruksanvisning | tr Kullanım bilgidir |

MADE IN GERMANY

BDKEY® Set

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Stand der Information 2016-12

BDKEY® Set

(en) Instructions for use, please read carefully!

Rx Only


Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Intended use

Besides functional impression-taking, the BDKEY® Set allows both the definition of the occlusal plane and the centre of face, as well as an aesthetic try-in combined with a functional check for fabricating full dentures. Furthermore, the BDKEY® Set serves to transfer analog data to the digital area of the Baltic Denture System.

BDKEY® Set consists of

- 3 pairs Upper- / LowerKEY PEEK in the sizes S, M and L
- 2 x BDKEY® Plane
- 2 x BDKEY® Fin
- 2 x BDKEY® Connect
- 10 x BDKEY® Lock

Individually adaptable BDKEY®s (made of PMMA) are available separately for patient-oriented work

BDKEY® PMMA (for single use and recommended for individual adaptation to the patient situation)

Indications

- **Upper- / LowerKEY (made of PEEK)** serve to determine the patient-individual size of the PMMA BDKEY® to be used.
- **Upper- / LowerKEY PMMA / BDKEY® Lock** serve for functional impression-taking, for checking the aesthetics of the rows of teeth (try-in) and for determining the position of the occlusal plane, as well as vertical and horizontal interlocking of UpperKEY / LowerKEY using BDKEY® Lock (determination of maxillomandibular relationship). BDKEY® (made of PMMA) is characterised by good and easy processing with rotary instruments and can be adapted to the individual patient situation during use by grinding at any time.

Note:

Sufficient suction of the impression, especially in the upper jaw is a prerequisite for the use of BDKEY® Set.

- **BDKEY® Plane / BDKEY® Fin** serve to align the UpperKEY, taking the occlusal plane and the patient's centre of face into consideration, as well as other aesthetic aspects.

Indication

- Edentulism of both jaws with a eugnathic jaw situation (Angle Class I)

Contraindication

- Partially edentulous jaw situation
- Dysgnathia (Angle Class II and III)
- Hypersensitivity to one of the components.

Composition

- UpperKEY / LowerKEY
medical polyetheretherketone (PEEK) pigmented with titanium dioxide
- ^{BD}KEY® Lock
medical copolymer based on polymethylmethacrylate (PMMA)
- ^{BD}KEY® Plane / ^{BD}KEY® Fin / ^{BD}KEY® Connect
medical polyetheretherketone (PEEK), carbon fibre-reinforced
- ^{BD}KEY® PMMA
medical copolymer based on polymethylmethacrylate (PMMA)

Use is only admissible by trained technical personnel.

Prior to use

- By trying-in and adapting the UpperKEY and LowerKEY (made of PEEK), the patient-specific applicability of the Baltic Denture System has to be checked beforehand and the required size of the medical copolymer based on polymethylmethacrylate (PMMA) determined.
- The occlusion and connector surfaces of UpperKEY and LowerKEY have to be checked prior to using adhesives (impression compound etc.).
- Instructions for sterilisation - see chap. 'Cleaning, visual inspection and sterilisation'.
- The flatness of ^{BD}KEY® Plane has to be ensured prior to use.
- Only use ^{BD}KEY® Set with the patient in an upright posture.

Use of the ^{BD}KEY® Set is divided into the following work steps

- Step 1: Assembly of ^{BD}KEY® Plane / ^{BD}KEY® Fin / ^{BD}KEY® Connect
- Step 2: Try-in and determination of the ^{BD}KEY® size

Note:


UpperKEY and LowerKEY of different sizes are not compatible. It is recommended to use the ^{BD}KEY® (made of PEEK) to determine the size, as it is reusable following reconditioning.

- Step 3: Individualisation and use of the UpperKEY PMMA
- Step 4: Individualisation and use of the LowerKEY PMMA
- Step 5: Investigation of function

Step 1

Assembly of ^{BD}KEY® Plane / ^{BD}KEY® Fin / ^{BD}KEY® Connect

- ^{BD}KEY® Fin is pushed centrally on the ^{BD}KEY® Plane into the groove provided, see sketch.
- ^{BD}KEY® Connect is attached to the ^{BD}KEY® Plane according to the sketch.

Step 2

Try-in and determination of the ^{BD}KEY® size

- Size selection of UpperKEY PMMA takes into consideration the patient's anatomical conditions from a try-in of ^{BD}KEY® PEEK. (Measurement of the old denture / alveolar ridge may be helpful.)
- Anatomical reference points then have to be marked (vertical height).

Step 3

Individualisation and use of the UpperKEY PMMA

- The size of the UpperKEY PMMA may be adapted with rotary instruments (e.g. recessing lip frenulum, position of the palatal vibrating line).
- The use of suitable primers is recommended to ensure sufficient adhesion of the impression material to the underside of the UpperKEY PMMA.
- To align the UpperKEY PMMA, it is recommended to place basal stops made of ^{BD}Impress® or a suitable thermoplastic material.
- While curing the stops, the rows of teeth of the UpperKEY PMMA have to be aligned with the aid of ^{BD}KEY® Plane and ^{BD}KEY® Fin according to functional aspects (Camper's plane, bipupillary line) as well as aesthetic considerations (centre of face, length of the anterior teeth).
- Build-up of the functional ridge is achieved in one or more work steps with ^{BD}Impress® or a comparable thermoplastic material. For functional impression-taking, the prepared UpperKEY PMMA is filled with suitable impression material and after repositioning is shaped under functional movement.
- The UpperKEY PMMA has to be removed and prepared for subsequent correction impression-taking.
- A suitable bonding agent has to be used prior to each use of impression materials.

Step 4

Individualisation and use of the LowerKEY PMMA

- UpperKEY PMMA and LowerKEY PMMA are connected using ^{BD}KEY® Lock.
- The interlocked UpperKEY and LowerKEY are then placed in the patient's mouth to check the spatial conditions between the alveolar ridge and the basal surface of LowerKEY PMMA. Any early contacts of the basal surface have to be eliminated by grinding.
- After placing three stops in the basal area of the LowerKEY PMMA using ^{BD}Impress® or a comparable thermoplastic material, the patient moves their lower jaw into its familiar, centric terminal occlusion. Determination of the relation between the upper and lower jaws in terms of the vertical and horizontal position (intercondylar axis) is thus achieved.
- The vestibular functional ridge is built up with ^{BD}Impress® or a comparable thermoplastic material and is shaped under functional movement in the mouth. Following separation of the ^{BD}KEY® PMMA, the lingual functional ridge is built up and shaped in the mouth while moving the tongue.
- It is useful to undertake an intermediate check of the maxillomandibular relationship prior to fine impression-taking. In case of significant deviations, the previous steps to determine the relation between the vertical and horizontal positions of the upper and lower jaws have to be repeated and corrected.

- The UpperKEY PMMA / LowerKEY PMMA should remain interlocked and the entire surface filled with low viscosity impression material for fine impression-taking. Impression-taking of the functional ridge is performed while undertaking facial movements.
- A suitable bonding agent has to be used prior to each use of impression materials.

Step 5

Investigation of function

- After unlocking UpperKEY and LowerKEY by removing ^{BD}KEY® Lock, speech checks and functional movements in the mouth can be undertaken. It is also necessary to check centricity and aesthetics.
- UpperKEY and LowerKEY can be removed from the mouth.

Storage

Protect against direct sunlight.

Special notes

- Improper use of ^{BD}KEY® Set and ^{BD}KEY® PMMA may lead to the patient's aesthetic and functional impairments.
- Only use ^{BD}KEY® Set and ^{BD}KEY® PMMA with the patient in an upright posture.
- Both Upper- and LowerKEY PEEK, as well as PEEK ^{BD}KEY® Plane and ^{BD}KEY® Fin (made of PEEK), are sterilisable and therefore reusable.
- Reuse of the single-use items ^{BD}KEY® Lock and ^{BD}KEY® PMMA is inadmissible.
 - This can lead to inaccuracy of fit with ^{BD}KEY® Lock.
 - A contaminated ^{BD}KEY® Lock and ^{BD}KEY® PMMA may constitute a health hazard.
 - ^{BD}KEY® Lock and ^{BD}KEY® PMMA are non-sterilisable.
- The use of non-authorised system components in conjunction with the Baltic Denture System is inadmissible.
- The list of approved system components is available at www.baltic-denture-system.de in the **Service** section.
- Do not use products after the specified expiration date.

Reconditioning (^{BD}KEY® PEEK, ^{BD}KEY® Plane / -Fin / -Connect)

With the exception of ^{BD}KEY® Lock and ^{BD}KEY® PMMA (single use items), all components are visually inspected and sterilised prior to use as follows.

Cleaning and visual inspection

- Rinse off soiling under running water; take special care given the presence of cavities.
- Cleaning can be carried out with a soft brush.
- No additional cleaning agents should be used.
- After cleaning, the products should be visually inspected for residues or damage.

Sterilisation

(Validated method in accordance with ISO 14937 and ISO 17665-1)

- Select a package that is suitable for the product and sterilisation procedure.
- The package must be large enough that the seal is not under tension.
- Transparent sterile packages with welded or bonded seams:
 - ▶ Sealed seams must be at least 8 mm wide
 - ▶ Allow a distance at least 10 mm between the sealing seam and the insertion
- Parameters for steam sterilisation:
 - ▶ Fractionated pre-vacuum (Type B)
 - ▶ Sterilisation temperature: 134 °C / 273,2 °F
 - ▶ Hold time: at least 3 min (full cycle)
 - ▶ Drying time: at least 10 min
 - ▶ A sterilisation indicator should be on every sterilisation package, or preferably in every package.
- The labelling of the sterilisation package should include the following:
 - ▶ Content
 - ▶ Sterilisation and "Use by" date
 - ▶ Responsible person (acc. to Quality Management)

Disinfection (BDKEY® PMMA und BDKEY® Lock)

On conclusion of the dental work steps, the BDKEY® PMMA should be disinfected with impression-taking disinfectants that contain no alcohol or aldehydes.

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de.

Disposal

The BDKEY® Set and BDKEY® PMMA components are water-insoluble, inactive, do not pose a threat to groundwater and may therefore be disposed of with normal household waste.

The product properties are based on compliance with the instructions for use.

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Merz Dental GmbH