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
M-ZR white HT/color HT/multicolor ST/multilayer HT+

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
Sintering Instructions

M-ZR white HT/color HT/multicolor ST/multilayer HT+

The listed firing curves are recommendations. Since the measurement of the actual temperature in each furnace can lead to different results, it may be necessary in individual cases to adjust the individual furnace parameters by means of a test firing. We are happy to advise you in this respect..

Schematic representation of the sintering processes

Green Body

The blank after pressing
The organic binder is still in the structure.

Blanc - White Body

The porosity has closed to the desired level. The grain boundaries are evenly distributed. The blank now has very good machinability and edge stability. The porosity achieved supports the infiltration of coloring liquids.

Intermediate Stage

The porosity closes and the structure loses volume.

Final Sintered

After sintering, the material receives its final density ($> 6.0 \text{ g / cm}^3$) and strength. The final temperature, rate of increase and holding times have a considerable influence on the strength and durability of the frames.

Dental frameworks and crowns made of zirconia should be fired in a furnace which is only used for these products. If restorations made of other ceramic materials or liquid ceramics are fired in the same furnace, this can impair the sintering behaviour or cause localised discoloration.

Please also note the detailed Merz Dental processing instructions for our respective products in our download center on our homepage: www.merz-dental.de

I. Recommendations prior to sintering

1. "Sinterdrops"

If you have the option to apply "sinterdrops" on the occlusal surface, this is definitely recommended.

Drops serve to stabilise frameworks during the sintering process.

At least four drops must be placed to ensure a safe support of frameworks.

- Single copings do not need any sintering drops.
- Anterior bridges with three units don't need drops. Posterior bridges should always be designed with drops.
- End units should always have drops.
- Pontics and molar copings must always have drops.
- Premolar and anterior copings can be left out, however at least every second unit should have a drop.
- Drops must be applied alternating in the oral and vestibular region.
- Drops should not be placed directly on the edge of the occlusal surface.

2. Cutting out / separating the work

(M-ZR white HT/color HT/multicolor ST)


Bridges up to 7 units:

Constructions that have a span of up to seven units can be separated directly from the blank. It is sufficient to use a turbine without water cooling; work should be carried out at a low speed and with sharp cutting tools. In order to prevent tension and thus premature damage to the work, each retaining bar should only be cut halfway through. Subsequently cut bars until predetermined breaking point before separating the construction. After the work has been carefully removed from the blank, the protruding retaining bars can be sanded down.



Bridges of 8 and more units:

In order to achieve optimum and distortion-free sinter results large spanned bridges of 8 and more units must not be completely removed from the blank. Depending on the curve of the bridge (arc), the oral holding bars must be preserved here.

- The holding bars should only remain attached to matching units:
Bars should either remain attached to copings or to pontics.
- If possible, the holding bars should remain to the end units.
If the end units are not the same type, then the bar should be attached to the next matching unit.
- Two holding bars should remain on the end units if the bridge is almost symmetrical or if it only extends over one quadrant.
If a bridge is constructed asymmetrically, one holding bar is retained on one unit and the third on a unit in the dental arch.

If you leave a solid piece of the blank, you should ideally thin it out.

situation	schematic drawing	regulation
12 units, symmetrical		2 bars, end position, oral
11 units, almost symmetrical		2 bars, end position, oral
10 units, asymmetrical		2 bars, end position, oral 1 bar within the arch, oral
9 units, asymmetrical		1 bar, end position, oral 1 bar at the coping next to the cantilever bridge, oral 1 bar within the arch, oral
8 units, asymmetrical		2 bars, end position, oral 1 bar within the arch, oral
unilateral		2 bars, end position, oral

3. Dry firing, if milling or sanding was carried out wet

If "wet" processing was carried out, a drying firing should be carried out in a (ceramic) furnace in order to remove the cooling and / or lubricating fluid from the porous structure.

After milling and before further processing, the restorations should be freed of milling dust. For this, the framework can be briefly steam cleaned or carefully cleaned under water. Dry firing should then be carried out to remove water and milling additives from the porous structure. Dry firing is required to achieve a homogeneous colour result when using coloring liquid.

We recommend: Drying of frameworks at 80 °C (for 30 min.) or 150 °C (for 10 min.)

4. Coloring with M-ZR Coloring-Liquids:

Do not color wet constructions, only completely dry frames should be colored.

Please adhere to the specified drying times.

5. Inspection of the milling work

- No material outbreaks
- No identifiable cracks
- No adherence of zirconium dioxide dust (removal using compressed air, brush and / or immersion in water with subsequent drying in a convection oven)
- No impurities
- No surface discoloring
- No glossy patches on the surface

If you discover one of the listed defects, the construction must not be used to manufacture dentures.

II. Sintering process

1. General information

In addition to using the correct sintering program, including final temperature, holding times and heating and cooling rates, the quality of the sintered restorations can depend on the following factors:

- the firing capacity and chamber volume of the furnace used
- the size and volume of individual constructions
- the mass in the firing chamber (sintering aids and number of constructions).

In order to achieve ideal results, a sintering program should be selected that heats up all the constructions and parts in the furnace uniformly. All bridges require a somewhat longer time to heat up uniformly and for temperature equalisation, due to the differing material thickness of (abutment) crowns and bridge elements. Local differences in sintering caused by heating too rapidly can lead to warping and the formation of cracks. This effect is particularly observed in the case of long-span and solid bridges. Where your furnace is loaded with numerous parts and additional elements such as hoods or similar, the thermal energy may, depending on the furnace type involved, prove inadequate to ensure flawless sintering of all the parts. Solid covers can absorb large quantities of energy which is then no longer available for dental work. Sintering without a cover can also improve the optical properties.

Therefore, we recommend to sinter the zirconium constructions without cover!

General rule:

A sintering program with slower rising rates is always an advantage for an optimal quality of the end product.

2. Final temperatures

All mentioned final temperatures are recommendations.

The actually reached temperature varies from sintering furnace to sintering furnace.

An opaque result after sintering can be an indication of a wrong furnace temperature.

M-ZR multicolor ST and M-ZR multilayer HT+ → Final temperature 1450 °C

M-ZR color HT and M-ZR white HT → Final temperature 1450 °C - 1530 °C

Important:

Do not place milled work in the furnace where the maximum residual temperature is greater than 70°C (temperature shock).

Never open the furnace door at a temperature exceeding 200°C.

(Possible damage to heating elements and the objects being fired - temperature shock).

3. Sintering Programs for all product versions

3.1 Standard Program:

Maximum temperature 1450°C (up to 1530°C for M-ZR white HT and M-ZR color HT) / Duration ~ 9 h.

- For fully anatomical crowns and bridges, caps and bridge constructions without sintering support.
- Normal furnace filling (without cover).
- Program line C0 L0 T008.A900 L9 T1800 T003.C1450 T7200 T010.C200 C0 L0 T2

	Temp. 1 [°C]	Temp. 2 [°C]	Heating rate [°C/h]	Heating rate [°C/min]	Hold time [min]	Time [min]
Heating	20	900	480	8	-	110
Hold	900	900	-	-	30	30
Heating	900	1450	200	3	-	183
Hold	1450	1450	-	-	120	120
Cooling*	1450	200	600	10	-	125
Total time:						568 min. 9.4 h

***Alternative:**

Uncontrolled cooling down by turning off the heating. Do not open the furnace before 200°C.

3.2 Long Program:

Maximum temperature 1450°C (up to 1530°C for M-ZR white HT and M-ZR color HT) / Duration ~ 16 h.

- For fully anatomical crowns, solid bridges and frames with sintering support, e.g. blank-segment.
- High furnace filling.
- Program line C0 L0 T008.A900 L9 T1800 T002.C1450 T7200 T010.C200 C0 L0 T2

	Temp. 1 [°C]	Temp. 2 [°C]	Heating rate [°C/h]	Heating rate [°C/min]	Hold time [min]	Time [min]
Heating	20	900	480	8	-	110
Hold	900	900	-	-	30	30
Heating	900	1450	100	2	-	275
Hold	1450	1450	-	-	120	120
Cooling*	1450	200	600	10	-	125
Total time:						660 min. 11 h

***Alternative:**

Uncontrolled cooling down by turning off the heating. Do not open the furnace before 200°C.

4. Speed-Sinter Program:

Dekema Austromat 664					
Furnace filling: Only sinter on one level, a maximum of 3 crowns per sintering process					
Program line	C0 L0 T060.A990 L9 T010.C1350 T015.C1450 T3000 T010.C1350 C990 C0 L0 T2				
	Temp. 1 [°C]	Temp. 2 [°C]	Heating rate [°C/min]	Hold time [min]	Furnace lift position
Heating	20	990	60	-	completely closed
	990	1350	10	-	-
	1350	1450	15	-	-
Hold	-	1450	-	50	-
Cooling	1450	1350	10	-	-
	1350	990	quickest possible	-	-
	-	990	-	-	completely open, heating off

Dekema Austromat 674					
Furnace filling: Only sinter on one level, a maximum of 6 crowns per sintering process					
Program line	C0 L0 T060.A990 L9 T010.C1350 T015.C1450 T4800 T010.C1350 C990 C0 L0				
	Temp. 1 [°C]	Temp. 2 [°C]	Heating rate [°C/min]	Hold time [min]	Furnace lift position
Heating	20	990	60	-	completely closed
	990	1350	10	-	-
	1350	1450	15	-	-
Hold	-	1450	-	80	-
Cooling	1450	1350	10	-	-
	1350	990	quickest possible	-	-
	-	990	-	-	completely open, heating off

Attention:

The specified speed-sinter cycles were validated in the Dekema Austromat 664/674 furnaces mentioned above. The use of the programs in other sintering furnaces can lead to different results due to different chamber volume, especially with regards to color and translucency of the constructions.

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