PEEK BioSolution®

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(en) Instructions for use, please read carefully!

Indication

- Complete loss of individual teeth
- Abnormal jaw function
- Bruxism

Intended purpose

PEEK BioSolution®, to produce temporary and permanent dental restorations

- Hybrid abutments
- Gingiva formers
- Crowns and bridges
- Transversal connectors
- Sublingual brackets
- Retaining elements: Telescopic crowns, Attachments, Retaining clips, Bar constructions
- Frameworks for clasp dentures
- Denture bases
- Occlusal splints

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Notes for the user

Construction

- The minimum wall thickness for crowns and bridges should not be less than 0.7 mm.
- No more than two pontics may be planned between two posts.
- The connector cross-section should be maximised as far as possible in order to ensure maximum stability.
- Given two pontics in the posterior region, additional lingual or palatal reinforcement is recommended.
- Only parallel-walled bases preferably with circular support, are to be used in the manufacture of hybrid abutments.

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Veneering of PEEK frameworks

- PEEK BioSolution® can be individualised using common veneer and composite materials, and following conditioning can be veneered with the artConcept® PLUS system in conjunction with artVeneer® / artVeneer® life shells.
- Ideally a veneer resin is used with the same or similar low elastic modulus as the framework material (5.100 MPa). This is necessary to avoid the formation of gaps or detachment of very brittle resins (e.g. highly filled composite) due to the intrinsic mobility of the construction.
- The following veneer systems have been successfully tested with PEEK BioSolution®:
 - ▶ artConcept® PLUS (Merz Dental)
 - ▶ in:Joy (Dentsply DeguDent)*
 - ▶ signum® (Heraeus Kulzer)*
 - ▶ Solidex (Shofu)*

If other veneer systems are used, always ensure that especially the elastic modulus of the veneer composite lies in the same range as PEEK BioSolution® and a suitable primer is always used - please contact the manufacturer of the respective veneer system in this case.

Sterilisation with heat / steam

(Validated method in accordance with ISO 14937 and ISO 17665-1)

When using PEEK BioSolution® as a hybrid abutment, cervical gingiva former or for other subgingival indications or in the vicinity of damaged tissue, the construction has to be sterilised prior to according to the following procedure.

- Rinse off soiling under running water; take special care given the presence of drill holes.
- Select a package that is suitable for the product and sterilisation procedure.
- The package must be large enough that the seal is not under tension.
- Transparent sterile packages with welded or bonded seams:
 - ▶ Sealed seams must be at least 8 mm wide
 - ▶ Allow a distance at least 10 mm between the sealing seam and the insertion

Parameters for steam sterilisation:

- Fractionated pre-vacuum (Type B)
- Sterilisation temperature: 134°C
- Hold time: at least 3 min (full cycle)
- Drying time: at least 10 min
- A sterilisation indicator should be on every sterilisation package or preferably in every package.
- The labelling of the sterilisation package should include the following:
 - ▶ Content
 - ▶ Sterilisation and "Use by" date
 - ▶ Responsible person (acc. to Quality Management)

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- Protect from sunlight.

Notes

- Do not use the material after the specified expiration date.
- To ensure traceability and assignment of the product data as well as the LOT numbers for the patient at any time, it is recommended to store the milling blank in the original packaging during two milling operations.

- Responsibility has to be taken that the treatment (sterilisation) performed with the actual equipment used, materials and personnel achieves the required results. This normally requires routine monitoring of the validated mechanical procedure. In addition, any deviation from the procedure stated in these instructions for use should be carefully evaluated by the user for its effectiveness and possible adverse effects.
- Changes to the construction through use, e.g. due to fissures, cracks, clouding or significant signs of abrasion, necessitate a repair or refabrication. Monitoring should be carried out with regular assessment of the denture at the follow-up presentation in the dental surgery.
- All serious incidents occurring with users and / or patients in connection with this product must be reported to Merz Dental, as well as to the German Federal Institute for Drugs and Medical Devices (BfArM), or the respective competent authority of the member state in which the user and / or patient is resident.
- For more information please visit www.merz-dental.de/service

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For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de.

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- * The company names and names of third-party products are trademarks or registered trademarks of the respective manufacturer. The use of these names without reference to a (registered and / or protected) trademark or band mark etc. is not to be regarded as violation of the protective rights of these names and not as damage towards those companies, which possess these rights.

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PEEK BioSolution®

Milling blank made of inorganically reinforced PEEK

usa Instructions for use, please read carefully!

Rx Only!

(Caution: United States Law restricts this device to sale by or on the order of a licensed dental practitioner)

Intended use

PEEK BioSolution® is intended to be used for the fabrication of permanent dental restorations using CAD/CAM techniques. These are:

- Fully anatomical crowns and bridges (max. 2 pontics and min. 13 mm² connector cross-section)
- Crown copings and bridge substructures for composite veneering (max. 2 pontics and min. 13 mm² connector cross-section)
- Telescopic primary and secondary crowns and frameworks
- Secondary bar structures on primary bars made of titanium alloy, CoCr alloys, gold alloy, zirconium dioxide

Classification according to MDD 93/42/EEC and MDR (EU) 2017/745

Class IIb medical device for long-term contact with intact mucous membranes or temporary contact with injured mucous membranes.

Patient target groups

All genders and age groups

Intended users

Dental professionals

Contraindication

- Hypersensitivity to one of the components

Composition for white and A2/B2

Polyether ether ketone (PEEK)

< 20 % titanium dioxide (TiO₂)

< 1 % TiO₂ based pigment (only for shade A2/B2)

Symbols

Reference that the product is a medical device

Physical data (white & A2/B2)

Flexural strength in accordance with EN ISO 20795-1	170 MPa
Elastic modulus / flexural modulus in accordance with EN ISO 20795-1	5,100 MPa
Flexural strength in accordance with EN ISO 10477	200 MPa
Charpy notched impact strength 1eA in accordance with EN ISO 179	4.05 kJ/m ²
Charpy impact strength 1eU in accordance with EN ISO 179	> 99,81 ± 0 kJ/m ²
no material failure	
Vickers hardness acc. EN ISO 6507-1	32 HV 0.2
Ball indentation hardness EN ISO 2039-1 after 10 sec loading at 961 N	385 MPa
Total fracture work in accordance with EN ISO 20795-1	1,900 J/m ²
Standard testing for HI resins requires at least 900 J/m ²	
Fracture toughness in accordance with EN ISO 20795-1	4.7 MPa x √m
Standard testing for HI resins requires at least 1.9 MPa x √m	
Linear thermal expansion coefficient	1.7 × 10 ⁻⁵ /K
Water absorption in accordance with ISO 62	0.4 %
Sliding friction coefficient	Al ₂ O ₃ μf = 0,126 ZrO ₂ μf = 0,125
Cytotoxicity test in accordance with EN ISO 10993-5	Max. value for biocompatibility (= 100 %)

Notes for the user

Construction

- The minimum wall thickness for crowns and bridges should not be less than 0.7 mm.
- No more than two pontics may be planned between two posts.
- The connector cross-section should be maximized as far as possible in order to ensure maximum stability.
- Given two pontics in the posterior region, additional lingual or palatal reinforcement is recommended.

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CAM processing

- In dry processing in the milling machine, care must be taken that milling chips are extracted or evacuated.
- In wet processing, the processor must ensure that the coolants (additives) used are completely removed from the surface of the workpiece.
- PEEK BioSolution® can generally be used with the common milling strategies available for thermoplastics and PMMA. In case of questions, please contact the manufacturer of your milling machine.
- PEEK BioSolution® can generally be used with the common burs available for thermoplastics and PMMA. In case of questions, please contact the manufacturer of your burs.
- It is important, particularly for precise fitting parts (e.g. both between primary and secondary parts, as well as matrices) that sharp, preferably new burs are used. The use of single edge tools is recommended, whereby, for producing high precision fits, they should be replaced after 3 to 4 units at the latest.

- In comparison with other thermoplastics, PEEK BioSolution® has a very high softening temperature of ~ 345°C. It should nevertheless be ensured that there is no excessive heat transfer to the material during the milling process, e.g. caused by too high speeds and/or too high a feed rate.
- Punctual, intense overheating can discolour the material or lead to local microstructural damage. Slight deformations may otherwise occur, especially with very thin walled parts (e.g. friction caps of 0.4 mm - 0.5 mm).
- In the case of fitting parts, such as bar attachments, it is important in the final finishing operation to work with very small infeed steps (e.g. 0.05 mm) to obtain the smoothest possible surface.
- In the case of finely tapered parts that could evade the milling pressure, the first step recommended is to undertake finishing on the fitting side, which is not intended to be further processed (e.g. inner surfaces of braces or fraction caps). The outer side is processed in the second step.

Manual trimming

- Always:
Use sharp burs. The use of blunt or heavily worn burs is contraindicated for removing PEEK without damaging the material.
- Work with low contact pressure and moderate speeds (10,000 - 20,000 rpm).
- It is imperative to avoid initial damage to the material, e.g. with cutting disks while finishing saddle margins or interdental connectors etc. Otherwise the fracture and impact toughness of the material is significantly reduced due to microscopic cracks and notches.
- Fine processing is performed with conventional, cross-cut carbide burs.
- In case of thin construction parts, e.g. friction caps, vaporisation should be avoided if possible or special care must be taken.
- Prepolishing with low contact pressure using diamond or corundum-reinforced rubber polishers.
- High-gloss polish with common polishing pastes.

Surface conditioning prior to bonding / veneering / cementing

- Note PEEK BioSolution® is a highly compressed material that is largely resistant to chemical influences. The surface is highly sealed and smooth, the microstructure compact and there are no chemically reactive groups on the surface. This impedes the formation of chemical compounds with other polymers. Optimal surface conditioning prior to bonding / veneering / cementing is therefore extremely important and indispensable.
- Firstly the grease-free surface is evenly blasted with 110 µm aluminium oxide under a pressure of 2 to 4 bar. The blasting agent Rocatec™-Pre (3M ESPE)*, for example, has proven to be well suited for this purpose.
- Bonding systems have to be used. One-off blasting is not sufficient for optimal bonding and / or to avoiding cracking.

The following materials have become established for conditioning the surface roughened in this way:

- Rocatec™ Plus mit ESPE™ Sil (3M ESPE)*
- Z-PRIME™ plus (BISCO)*
- visio.link (bredent)*
- CLEARFIL™ CERAMIC PRIMER (Kuraray)*
- In the case of especially fine structures, Rocatec™ Soft 30 µm and ESPE™ Sil (3M ESPE)* can be used directly without pre-blasting. However, the bonding strength drops by around 25% as compared with Rocatec™ Plus (3M ESPE)*

Veneering of PEEK frameworks

- PEEK BioSolution® can be individualised using common veneer and composite materials.
- Ideally a veneer resin is used with the same or similar low elastic modulus as the framework material (5.100 MPa). This is necessary to avoid the formation of gaps or detachment of very brittle resins (e.g. highly filled composite) due to the intrinsic mobility of the construction.

Please contact the manufacturer of the respective veneer system to check for compatibility with PEEK.

Fixation in the mouth

- The conditioned construction can be temporarily fixed with eugenol-free fixing material.
- The conditioned construction can be permanently fixed with adhesive fixing material, e.g. composite cements

The following materials are suitable:

- ▶ Panavia™ F 2.0 (Kuraray)*

Please contact the manufacturer of the respective cement system to check for compatibility with PEEK.

Storage

- Protect from sunlight.

Special notes

- Do not use the material after the specified expiration date.
- To ensure traceability and assignment of the product data as well as the LOT numbers for the patient at any time, it is recommended to store the milling blank in the original packaging during two milling operations.
- Responsibility has to be taken that the treatment (sterilisation) performed with the actual equipment used, materials and personnel achieves the required results. This normally requires routine monitoring of the validated mechanical procedure. In addition, any deviation from the procedure stated in these instructions for use should be carefully evaluated by the user for its effectiveness and possible adverse effects.
- Responsibility has to be taken that the treatment (sterilisation) performed with the actual equipment used, materials and personnel achieves the required results. This normally requires routine monitoring of the validated mechanical procedure. In addition, any deviation from the procedure stated in these instructions for use should be carefully evaluated by the user for its effectiveness and possible adverse effects.
- Changes to the construction through use, e.g. due to fissures, cracks, clouding or significant signs of abrasion, necessitate a repair or refabrication. Monitoring should be carried out with regular assessment of the denture at the follow-up presentation in the dental surgery.
- All serious incidents occurring with users and / or patients in connection with this product must be reported to Merz Dental, as well as to the German Federal Institute for Drugs and Medical Devices (BfArM), or the respective competent authority of the member state in which the user and / or patient is resident.
- For more information please visit www.merz-dental.de/service

Risk phrases

- Processing PEEK BioSolution® may generate dust that can irritate the eyes, skin and airways. Therefore, always ensure that the extraction system on your milling machine and at the workplace is working properly for individual post-processing and that you are using personal protective equipment.
- For intraoral applications, small constructions must be secured against aspiration and swallowing (especially during insertion and removal).
- The use of a rubber dam is recommended with hybrid abutments, gingiva formers and single crowns, as well as constructions of a similar size.

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de. s.

Disposal

- Residues from processing PEEK BioSolution® have to be disposed of as residual waste.

The product properties are based on compliance with these instructions for use.

* The company names and names of third-party products are trademarks or registered trademarks of the respective manufacturer. The use of these names without reference to a (registered and / or protected) trademark or band mark etc. is not to be regarded as violation of the protective rights of these names and not as damage towards those companies, which possess these rights.

The physical data may fluctuate in a range of $\pm 5\%$ around the mean values specified.

The users themselves are responsible for the use of the product. Merz Dental assumes no liability and / or guarantee for the use of components that are not part of the system and / or for components that have not been tested and for incorrect results, as the manufacturer has no influence on processing. Any claims for damages that may still arise relate exclusively to the value of our products.

Date of information 2021-08

Merz Dental GmbH

Kieferweg 1 24321 Lütjenburg, Germany

Tel + 49 (0) 4381 / 403-0

Fax + 49 (0) 4381 / 403-403

www.merz-dental.de

EN ISO 13485

CE 0482

PEEK BioSolution®

Milling blank made of inorganically reinforced PEEK

(can) Instructions for use, please read carefully!

Intended use

PEEK BioSolution® is intended to be used for the fabrication of temporary (6 months) dental restorations using CAD/CAM techniques. These are:

- Fully anatomical crowns and bridges
- Crown copings and bridge substructures for composite veneering
- Telescopic primary and secondary crowns and frameworks
- Secondary bar structures on primary bars made of titanium alloy, CoCr alloys, gold alloy, zirconium dioxide

Contraindication

- Hypersensitivity to one of the components

Composition

	White	A2/B2	GUM
• Polyether ether ketone (PEEK)	~ 80%	~ 80%	-
• Titanium dioxide (TiO ₂)	~ 20%	~ 20%	~ 5%
• TiO ₂ based pigment	-	< 1%	-
• Iron oxide, red	-	< 1%	-

Symbols

MD Reference that the product is a medical device

Physical data

	White & A2/B2	GUM
Flexural strength in accordance with EN ISO 20795-1	170 MPa	170 MPa
Elastic modulus / flexural modulus in accordance with EN ISO 20795-1	5,100 MPa	4,500 MPa
Flexural strength in accordance with EN ISO 10477	200 MPa	220 MPa
Charpy notched impact strength 1eA in accordance with EN ISO 179	4.1 kJ/m ²	3.7 kJ/m ²
Charpy impact strength 1eU in accordance with EN ISO 179 no material failure	> 100 kJ/m ²	> 100 kJ/m ²
Vickers hardness acc. EN ISO 6507-1	32 HV 0.2	27 HV 0.2
Ball indentation hardness EN ISO 2039-1 after 10 sec loading at 961 N	385 MPa	350 MPa
Total fracture work in accordance with EN ISO 20795-1 Standard testing for HI resins requires at least 900 J/m ²	1,900 J/m ²	8,500 J/m ²
Fracture toughness in accordance with EN ISO 20795-1 Standard testing for HI resins requires at least 1.9 MPa x √m	4.7 MPa x √m	6.7 MPa x √m
Linear thermal expansion coefficient	1.7 × 10 ⁻⁵ /K	1.7 × 10 ⁻⁵ /K
Water absorption in accordance with ISO 62	0.4 %	0.4%
Sliding friction coefficient	Al ₂ O ₃ μf = 0,13 ZrO ₂ μf = 0,13	Al ₂ O ₃ μf = 0,09 ZrO ₂ μf = 0,10
Cytotoxicity test in accordance with EN ISO 10993-5	Max. value for biocompatibility (= 100 %)	

Notes for the user

Construction

- The minimum wall thickness for crowns and bridges should not be less than 0.7 mm.
- No more than two pontics may be planned between two posts.
- The connector cross-section should be maximised as far as possible in order to ensure maximum stability.
- Given two pontics in the posterior region, additional lingual or palatal reinforcement is recommended.

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CAM processing

- In dry processing in the milling machine, care must be taken that milling chips are extracted or evacuated.
- In wet processing, the processor must ensure that the coolants (additives) used are completely removed from the surface of the workpiece.
- PEEK BioSolution® can generally be used with the common milling strategies available for thermoplastics and PMMA. In case of questions, please contact the manufacturer of your milling machine.
- PEEK BioSolution® can generally be used with the common burs available for thermoplastics and PMMA. In case of questions, please contact the manufacturer of your burs.
- It is important, particularly for precise fitting parts (e.g. both between primary and secondary parts, as well as matrices) that sharp, preferably new burs are used. The use of single edge tools is recommended, whereby, for producing high precision fits, they should be replaced after 3 to 4 units at the latest.
- In comparison with other thermoplastics, PEEK BioSolution® has a very high softening temperature of ~ 345°C. It should nevertheless be ensured that there is no excessive heat transfer to the material during the milling process, e.g. caused by too high speeds and/or too high a feed rate.
- Punctual, intense overheating can discolour the material or lead to local microstructural damage. Slight deformations may otherwise occur, especially with very thin walled parts (e.g. friction caps of 0.4 mm - 0.5 mm).
- In the case of fitting parts, such as bar attachments, it is important in the final finishing operation to work with very small infeed steps (e.g. 0.05 mm) to obtain the smoothest possible surface.
- In the case of finely tapered parts that could evade the milling pressure, the first step recommended is to undertake finishing on the fitting side, which is not intended to be further processed (e.g. inner surfaces of braces or fraction caps). The outer side is processed in the second step.

Manual trimming

- Always:
Use sharp burs. The use of blunt or heavily worn burs is contraindicated for removing PEEK without damaging the material.
- Work with low contact pressure and moderate speeds (10,000 - 20,000 rpm).
- It is imperative to avoid initial damage to the material, e.g. with cutting disks while finishing saddle margins or interdental connectors etc. Otherwise the fracture and impact toughness of the material is significantly reduced due to microscopic cracks and notches.
- Fine processing is performed with conventional, cross-cut carbide burs.
- In case of thin construction parts, e.g. friction caps, vaporisation should be avoided if possible or special care must be taken.
- Prepolishing with low contact pressure using diamond or corundum-reinforced rubber polishers.
- High-gloss polish with common polishing pastes.

Surface conditioning prior to bonding / veneering / cementing

- Note PEEK BioSolution® is a highly compressed material that is largely resistant to chemical influences. The surface is highly sealed and smooth, the microstructure compact and there are no chemically reactive groups on the surface. This impedes the formation of chemical compounds with other polymers. Optimal surface conditioning prior to bonding / veneering / cementing is therefore extremely important and indispensable.
- The grease-free surface is evenly blasted with 110 µm aluminum oxide under a pressure of 2 to 4 bar.
- Bonding systems have to be used.
One-off blasting is not sufficient for optimal bonding and / or to avoiding cracking.

The following materials have become established for conditioning the surface roughened in this way:

- Rocatec™ Plus mit ESPE™ Sil (3M ESPE)*
- Z-PRIME™ plus (BISCO)*
- visio.link (bredent)*
- CLEARFIL™ CERAMIC PRIMER (Kuraray)*
- In the case of especially fine structures, Rocatec™ Soft 30 µm and ESPE™ Sil (3M ESPE)* can be used directly without pre-blasting.
However, the bonding strength drops by around 25 % as compared with Rocatec™ Plus (3M ESPE)*.

Veneering of PEEK frameworks

- PEEK BioSolution® can be individualised using common veneer and composite materials.
- Ideally a veneer resin is used with the same or similar low elastic modulus as the framework material (5.100 MPa). This is necessary to avoid the formation of gaps or detachment of very brittle resins (e.g. highly filled composite) due to the intrinsic mobility of the construction.

Please contact the manufacturer of the respective veneer system to check for compatibility with PEEK.

Fixation in the mouth

- The conditioned construction can be temporarily fixed with eugenol-free fixing material.

Wear period

- Can be worn up to 6.0 months

Please contact the manufacturer of the respective cement system to check for compatibility with PEEK.

Storage

- Protect from sunlight.

Special notes

- Do not use the material after the specified expiration date.
- Record the LOT number with every process that requires identification of the material.
- Responsibility has to be taken that the treatment (sterilisation) performed with the actual equipment used, materials and personnel achieves the required results. This normally requires routine monitoring of the validated mechanical procedure. In addition, any deviation from the procedure stated in these instructions for use should be carefully evaluated by the user for its effectiveness and possible adverse effects.
- Changes to the construction through use, e.g. due to fissures, cracks, clouding or significant signs of abrasion, necessitate a repair or refabrication. Monitoring should be carried out with regular assessment of the denture at the follow-up presentation in the dental surgery.
- All serious incidents occurring with users and / or patients in connection with this product must be reported to Merz Dental, as well as to the German Federal Institute for Drugs and Medical Devices (BfArM), or the respective competent authority of the member state in which the user and / or patient is resident.
- For more information please visit www.merz-dental.de/service

Risk phrases

- Processing PEEK BioSolution® may generate dust that can irritate the eyes, skin and airways. Therefore, always ensure that the extraction system on your milling machine and at the workplace is working properly for individual post-processing and that you are using personal protective equipment.
- For intraoral applications, small constructions must be secured against aspiration and swallowing (especially during insertion and removal).

Disposal

- PEEK BioSolution® is water-insoluble, chemically and biologically inactive, does not pose a threat to groundwater and may therefore be disposed of with normal household waste.

The product properties are based on compliance with these instructions for use.

* The company names and names of third-party products are trademarks or registered trademarks of the respective manufacturer. The use of these names without reference to a (registered and / or protected) trademark or band mark etc. is not to be regarded as violation of the protective rights of these names and not as damage towards those companies, which possess these rights.

The physical data may fluctuate in a range of $\pm 5\%$ around the mean values specified.

The users themselves are responsible for the use of the product. Merz Dental assumes no liability and / or guarantee for the use of components that are not part of the system and / or for components that have not been tested and for incorrect results, as the manufacturer has no influence on processing. Any claims for damages that may still arise relate exclusively to the value of our products.

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