PEEK *NATURE*

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- | | | | |
|----|----------------------|----|---------------------|
| en | Instructions for use | sv | Bruksanvisning |
| fr | Notice d'utilisation | da | Brugsanvisningen |
| it | Istruzioni per l'uso | no | Brugsanvisningen |
| es | Instrucciones de uso | hu | Használati utasítás |
| pl | Instrukcja użycia | el | Οδηγίες χρήσης |
| nl | Gebruiksaanwijzing | tr | Kullanım kılavuzu |

PEEK NATURE

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Stand der Information 2018-11

PEEK NATURE

Milling blank made of unfilled polyetheretherketone

(en) Instructions for use, please read carefully!

Indication

- Total or partial tooth loss

Intended use

PEEK Nature (EU medical device Class IIa) is suitable for the temporary or permanent restoration of chewing function by replacing lost tooth substance using

- Single crowns / crown frames (caps)
- Bridges / bridge frames
- Transversal connectors
- Sublingual brackets
- Denture bases
- Retaining elements
 - Telescopic crowns
 - Attachments
 - Retaining clips
 - Bar constructions

Contraindication

- Hypersensitivity to one of the components

Composition

- 100 % polyetheretherketone (PEEK)

Physical data

Flexural strength in accordance with EN ISO 20795-1	164 MPa
Elastic modulus / flexural modulus in accordance with EN ISO 20795-1	4,100 MPa
Flexural strength in accordance with EN ISO 10477	182 MPa
Charpy notched impact strength 1eA in accordance with EN ISO 179	4.2 kJ/m ²
Charpy impact strength 1eU in accordance with EN ISO 179 no material failure	>500 kJ/m ²
Vickers hardness acc. EN ISO 6507-1	23 HV 0.2
Ball indentation hardness EN ISO 2039-1 after 10 sec loading at 961 N	230 MPa
Total fracture work in accordance with EN ISO 20795-1 Standard testing for HI resins requires at least 900 J/m ²	14,500 J/m ²
Fracture toughness in accordance with EN ISO 20795-1 Standard testing for HI resins requires at least 1.9 MPa x \sqrt{m}	7.00 MPa x \sqrt{m}
Linear thermal expansion coefficient	1.9 MPa x \sqrt{m}
Water absorption in accordance with ISO 62	0.4 %

Notes for the user

Construction

- The minimum wall thickness for crowns and bridges should not be less than 0.7 mm.
- No more than two pontics may be planned between two posts.
- The connector cross-section should be maximised as far as possible in order to ensure maximum stability.
- Given two pontics in the posterior region, additional lingual or palatal reinforcement is recommended.
- Only parallel-walled bases preferably with circular support, are to be used in the manufacture of hybrid abutments.

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CAM processing

- In dry processing in the milling machine, care must be taken that milling chips are extracted or evacuated.
- In wet processing, the processor must ensure that the coolants (additives) used are completely removed from the surface of the workpiece.
- PEEK *Nature* can generally be used with the common milling strategies available for thermoplastics and PMMA. In case of questions, please contact the manufacturer of your milling machine and/or your software provider.
- PEEK *Nature* can generally be used with the common burs available for thermoplastics and PMMA. In case of questions, please contact the manufacturer of your burs.
- It is important, particularly for precise fitting parts (e.g. both between primary and secondary parts, as well as matrices) that sharp, preferably new burs are used. The use of single edge tools is recommended, whereby, for producing high precision fits, they should be replaced after 3 to 4 units at the latest.
- In comparison with other thermoplastics, PEEK *Nature* has a very high softening temperature of ~ 345°C. It should nevertheless be ensured that there is no excessive heat transfer to the material during the milling process, e.g. caused by too high speeds and/or too high a feed rate.
- Punctual, intense overheating can discolour the material or lead to local microstructural damage. Slight deformations may otherwise occur, especially with very thin walled parts (e.g. friction caps of 0.4 mm - 0.5 mm).
- In the case of fitting parts, such as bar attachments, it is important in the final finishing operation to work with very small infeed steps (e.g. 0.05 mm) to obtain the smoothest possible surface.
- In the case of finely tapered parts that could evade the milling pressure, the first step recommended is to undertake finishing on the fitting side, which is not intended to be further processed (e.g. inner surfaces of braces or fraction caps). The outer side is processed in the second step.

Manual trimming

- Please use sharp milling tools to be gentle on the material.
- Work with low contact pressure and moderate speeds (10,000 - 20,000 rpm).
- It is imperative to avoid initial damage to the material, e.g. with cutting disks while finishing saddle margins or interdental connectors etc. Otherwise the fracture and impact toughness of the material is significantly reduced due to microscopic cracks and notches.
- Fine processing is performed with conventional, cross-cut carbide burs.
- In case of thin construction parts, e.g. friction caps, vaporisation should be avoided if possible or special care must be taken.
- Perform pre-polishing with low contact pressure using diamond or corundum-reinforced rubber polishers as well as high-gloss polishing with common polishing pastes.

Surface conditioning prior to bonding / veneering / cementing

- Note PEEK *Nature* is a highly compressed material that is largely resistant to chemical influences. The surface is highly sealed and smooth, the microstructure compact and there are no chemically reactive groups on the surface. This impedes the formation of chemical compounds with other polymers. Conditioning of the surfaces and the use of a bonding system before bonding / veneering / cementing is therefore necessary.
- Firstly the grease-free surface is evenly blasted with 110 µm aluminium oxide under a pressure of 2 to 2.5 bar. The blasting agent Rocatec™-Pre (3M ESPE)*, for example, has proven to be well suited for this purpose.

The following materials have become established for conditioning the surface roughened in this way:

- Rocatec™ Plus with ESPE™ Sil and Visio™ Bond (3M ESPE)*
- visio.link (bredent)*

Veneering of PEEK frameworks

- PEEK *Nature* can be individualised using common veneer and composite materials, and following conditioning can be veneered with the artConcept® system in conjunction with artVeneer® veneer shells.
- The following veneer systems have been successfully tested with PEEK *Nature*:
 - ▶ artConcept® (Merz Dental)

When using other veneering systems, please observe the respective manufacturer's instructions.

Fixation in the mouth

- The conditioned construction can be temporarily fixed with eugenol-free fixing material.
- The conditioned construction can be permanently fixed with adhesive fixing material, e.g. composite cements

The following materials have been successfully tested with PEEK *Nature*:

- ▶ Panavia™ F 2.0 (Kuraray)*

Storage

- Protect from sunlight.

Notes

- Record the LOT number with every process that requires identification of the material.

Warnings

- Processing PEEK *Nature* may generate dust that can irritate the eyes, skin and airways. Therefore, always ensure that the extraction system on your milling machine and at the workplace is working properly for individual post-processing and that you are using personal protective equipment.
- For intraoral applications, small constructions must be secured against aspiration and swallowing (especially during insertion and removal).

Disposal

- Residues from processing PEEK *Nature* have to be disposed of as residual waste.

The product properties are based on compliance with these instructions for use.

* The company names and names of third-party products are trademarks or registered trademarks of the respective manufacturer. The use of these names without reference to a (registered and / or protected) trademark or band mark etc. is not to be regarded as violation of the protective rights of these names and not as damage towards those companies, which possess these rights.

The physical data may fluctuate in a range of $\pm 5\%$ around the mean values specified.

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Merz Dental GmbH