
artBloc® Temp

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
Stand der Information 2019-07

artBloc[®] Temp

Monochrome tooth-coloured block made of highly cross-linked interpenetrated Organic Modified Polymer-Network - OMP-N[®]

en Instructions for use, please read carefully!

Intended use

For functional restoration by replacement of lost tooth substance in permanent dentition in the form of temporary crowns and bridges in the anterior and posterior regions.

Indication

Loss of hard tooth substance and resulting aesthetic impairment.

Contraindication

Hypersensitivity to one of the components

Composition

- Polymethylmethacrylate (PMMA) and cross-linked polymers based on methacrylic acid esters
- Colourants
- Residual peroxide as dibenzoyl peroxide
- Methylmethacrylate (MMA) may be contained as residual monomer up to max. 1%.

Notes for the user

- Minimum wall thickness after polishing ≥ 1 mm
- In the construction of bridges with one pontic, the connection cross-section between the crown and pontic must not be less than 9 mm².
- In the construction of bridges spanning up to two pontics, the connection cross-section between crowns and pontics must not be less than 12 mm².
- By selecting the respective fixation system according to the indication, the user has to decide case-by-case whether the OMP-N[®] (PMMA) material ensures sufficient chewing stability for the restoration.

PHYSICAL PROPERTIES

| | |
|----------------------|------------------------------------|
| Flexural strength | 93 MPa according EN ISO 10477 |
| Module of elasticity | 2,680 MPa according EN ISO 20795-1 |

Veneering of PMMA frameworks

- Milled PMMA frameworks, crowns and bridges may be coated with conventional veneering materials. Materials based on MMA are preferable here.
- When veneering with composites the connection with the PMMA must always be made according to the respective manufacturer's specifications.

Trimming

- Cross-cut carbide burs suitable for resin or suitable cutting discs should be used to rework and separate frameworks, crowns and bridges from the blank.
- To prevent inaccurate fitting after fabricating the construction, avoid excessive heat development during trimming and polishing.

Restoration

Frameworks, crowns and bridges which have been fabricated from the artBloc® Temp can be restored and supplemented at any time with all standard MMA-based cold-curing polymers. Prior conditioning of the surface is necessary, preferably with a standard MMA-based bonding agent.

Fixation in the mouth

for the temporary restoration: for long-term temporary restorations:
- temporary cements - glass ionomer cements without bonding

Wear period

Can be worn up to 3 years

Storage

No special measures required.

Notes

- Changes to the construction through use, e.g. due to fissures, cracks or significant signs of abrasion, necessitate a repair or refabrication. Monitoring should be carried out with regular assessment of the denture at the follow-up presentation in the dental surgery.
- Record the LOT number with every process that requires identification of the material.
- To ensure traceability and assignment of the product data as well as the LOT numbers for the patient at any time, it is recommended to store the milling blank in the original packaging during two milling operations.

Warnings

Processing PMMA frameworks may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system on your milling machine and at the workplace is working properly for individual post-processing and that you are using personal protective equipment.

Disposal

Residues from processing artBloc® Temp have to be disposed of as residual waste.


For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de.

The product properties are based on compliance with and observation of these instructions for use.

Date of information 2019-04

artBloc[®] Temp

Monochrome tooth-coloured block made of highly cross-linked interpenetrated Organic Modified Polymer-Network - OMP-N[®]

usa Instructions for use, please read carefully!

Note: The artBloc temp is designed for milling equipment from imes-icore, Organical, Zirkozahn or Equivalent

Caution: Federal (USA) law restricts this device by or on the order of a licensed dental practitioner.

Intended use

For use as a milling blank in the fabrication of provisional crowns and bridges.

Indication

Loss of hard tooth substance and resulting aesthetic impairment.

Contraindication

Hypersensitivity to one of the components

Composition

- Polymethylmethacrylate (PMMA) and cross-linked polymers based on methacrylic acid esters
- Colourants
- Residual peroxide as dibenzoyl peroxide
- Methylmethacrylate (MMA) may be contained as residual monomer up to max. 1%.

Notes for the user

- Minimum wall thickness after polishing ≥ 1 mm
- In the construction of bridges with one pontic, the connection cross-section between the crown and pontic must not be less than 9 mm².
- In the construction of bridges spanning up to two pontics, the connection cross-section between crowns and pontics must not be less than 12 mm².
- By selecting the respective fixation system according to the indication, the user has to decide case-by-case whether the OMP-N[®] (PMMA) material ensures sufficient chewing stability for the restoration.

PHYSICAL PROPERTIES

| | |
|----------------------|------------------------------------|
| Flexural strength | 93 MPa according EN ISO 10477 |
| Module of elasticity | 2,680 MPa according EN ISO 20795-1 |

Veneering of PMMA substructures

- A milled PMMA substructure can be veneered with conventional appropriate materials. It is recommended to use material on a MMA basis.
- For composite veneers, the bond to the PMMA must always be done according to the instructions of the manufacturer.

Trimming

- Cross-cut carbide burs suitable for resin or suitable cutting discs should be used to rework and separate crowns and bridges from the blank.
- To prevent inaccurate fitting after fabricating the construction, avoid excessive heat development during trimming and polishing.

Restoration

Crowns and bridges which have been fabricated from the artBloc® Temp can be restored and supplemented at any time with all standard MMA-based cold-curing polymers. Prior conditioning of the surface is necessary, preferably with a standard MMA-based bonding agent.

Fixation in the mouth

for the temporary restoration: for long-term temporary restorations:
- temporary cements - glass ionomer cements without bonding

Wear period

Can be worn up to 3 years

Storage

No special measures required. See label for expiration date.

Notes

- Changes to the construction through use, e.g. due to fissures, cracks or significant signs of abrasion, necessitate a repair or refabrication. Monitoring should be carried out with regular assessment of the denture at the follow-up presentation in the dental surgery.
- Record the LOT number with every process that requires identification of the material.
- To ensure traceability and assignment of the product data as well as the LOT numbers for the patient at any time, it is recommended to store the milling blank in the original packaging during two milling operations.

Risk phrases

- Processing PMMA discs may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system on your milling machine and at the workplace is working properly for individual post-processing and that you are using personal protective equipment.

Disposal

Residues from processing artBloc® Temp have to be disposed of as residual waste.


For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de.

The product properties are based on compliance with and observation of these instructions for use.

Date of information 2019-04

artBloc[®] Temp

Monochrome tooth-coloured block made of highly cross-linked interpenetrated Organic Modified Polymer-Network - OMP-N[®]

can Instructions for use, please read carefully!

Purpose of use

artBloc[®] Temp is a Monochrome tooth-coloured block made of highly cross-linked interpenetrated OMP[®]-N (Organic Modified Polymer Network)

- acrylic crowns, bridges and framework for provisional restorations within the preprosthetic treatment of anterior and posterior teeth.
- In constructing bridges the maximum number of units is 6, number of pontics 2. Cantilevered pontics: none.

Composition

- Polymethylmethacrylate (PMMA) and cross-linked polymers based on methacrylic acid esters
- Colourants
- Residual peroxide as dibenzoyl peroxide
- Methylmethacrylate (MMA) may be contained as residual monomer up to max. 1%.

Finishing

- For refinishing and removing the framework from the blank, it is recommended to use helical-tooth cutters (hard metal) or suitable cutting discs.
- In order to avoid problems with the marginal fit after constructing the framework, it is recommended to avoid excessive heat during the preparation and polishing.

Repair

Frameworks, whole crowns and bridges manufactured from the artBloc temp can be repaired and added on to at any time with any customary cold-cure polymer on MMA basis.

Options for oral attachment

| | |
|----------------------------------|---|
| For the provisional restoration: | for long-term temporaries: |
| - temporary cements | - glass ionomer cements without bonding |

Storage

No special measures required.

Special instructions

Indicate the LOT No. for every procedure requiring an identification of the material.

Safety precautions

Dust is created during processing, which may lead to irritations of the eyes, skin and airways. Therefore, always pay attention that the exhaust system of your milling machine is working perfectly.

Further safety information can be obtained from the material safety data sheet.

The product properties are based on compliance with and observation of this instruction for use.

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Merz Dental GmbH