# artOpaque**PLUS**

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- |                               |                             |
|-------------------------------|-----------------------------|
| (en) Instructions for use     | (it) Naudojimo instrukcija  |
| (fr) Notice d'utilisation     | (lv) Lietošanas informācija |
| (it) Istruzioni per l'uso     | (et) Kasutusjuhend          |
| (es) Instrucciones de uso     | (hu) Használati utasítás    |
| (pt) Instruções de utilização | (cs) Návod k použití        |
| (pl) Instrukcja użycia        | (sl) Navodila za uporabo    |
| (nl) Gebruiksaanwijzing       | (sk) Návod na používanie    |
| (sv) Bruksanvisning           | (bg) Инструкции за употреба |
| (da) Brugsanvisningen         | (el) Οδηγίες χρήσης         |
| (no) Brugsanvisningen         | (tr) Kullanım kılavuzu      |
| (fi) Käyttöohje               |                             |

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# artOpaque<sup>PLUS</sup>

---

Tooth-shade paste opaquer, light-curing.

**(en) Instructions for use, please read carefully!**

## Indication

- Total or partial tooth loss
- Loss of hard tooth substance

## Intended use

Replacement of lost hard tooth substance by:

- Tooth-shaded fixture and for backing artVeneer<sup>®</sup>
- Veneering telescopic, attachment, bar and implant constructions with artVeneer<sup>®</sup>
- Tooth-shaded backing of trimmed acrylic teeth
- Fixing acrylic teeth to the model cast with a back protection plate or pontic
- Individualising acrylic teeth, artVeneer<sup>®</sup> and CAD/CAM restorations produced from M-PM<sup>®</sup> Disc and the artBloc<sup>®</sup> Temp
- Repairs

## Contraindication

- Hypersensitivity to one of the components

## Components of the artConcept<sup>®</sup> PLUS system

- artPrime PLUS - Metal bonding agent
- artPreOpaque PLUS - Priming opaquer
- **artOpaque PLUS - Opaquer for covering framework**
- artDentine PLUS/artEnamel PLUS Polymer - Tooth-shade PMMA acrylic
- artDentine Liquid NT PLUS (normal time) / artDentine Liquid LT PLUS (long time) - Monomer component for artDentine PLUS / artEnamel PLUS
- artConnect PLUS - Acrylic bonder
- artVeneer<sup>®</sup> - Tooth-shade PMMA anterior and posterior veneer

## Composition

Butandiol dimethacrylate diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide

## Processing

Apply artOpaque PLUS twice thinly using a brush with short, stiff bristles until the framework is completely covered with colour.

**Tip:** By applying the opaquer thinly, the complete polymerisation of the individual layers is ensured.

## Trimming

- Cross-cut carbide burs suitable for resin or suitable cutting discs should be used to rework.
- To prevent inaccurate fitting after polymerisation, avoid excessive heat development during trimming and polishing.

## Curing unit artOpaque PLUS

DC-XS / Uni-XS	Kulzer	90 sec
HeraFlash	Kulzer	90 sec
Labolight LV-II / III	GC	1 min
MPA 200	Dentsply	1 x HD
Solidilite EX / V	Shofu	1 min
Spektra 2000	Schütz Dental	3 min
Spektra 2000 fast	Schütz Dental	3 min
Spektra LED	Schütz Dental	30 sec
Spektramat	ivoclar vivadent	2 min
Triad II	Dentsply	5 min
LicuLite	Dentsply	1 x HD
Lumamat 100	ivoclar vivadent	Program 2
Targis Power	ivoclar vivadent	2 min

## Storage

Store in a dry and cool area at 10 - 25 °C / 50 - 77 °F without direct exposure to light.

## Disposal

The fully polymerized components of artConcept® PLUS are disposed of as residual waste. Disposal of non-polymerized product residues according to legal regulations. Do not allow entry into the sewage system or waterways.

## Notes

- Close all containers tightly after use.
- Do not use the material after the specified expiration date.
- To reduce the residual monomer content, the finished (sanded) but not polished polymerised workpiece should be stored in water at room temperature for about 48 hours before delivery.
- The residual monomer content after polymerisation of artDentine PLUS/artEnamel PLUS is < 1.8%.

## Danger and safety instructions (Contents <125ml)

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P350 IF ON SKIN: Gently wash with plenty of soap and water.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 Take off contaminated clothing and wash it before reuse.

---

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from [www.merz-dental.de](http://www.merz-dental.de).

---

**The product properties are based on compliance with and observation of these instructions for use.**

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