# artPreOpaque**PLUS**

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- |                                    |                                       |
|------------------------------------|---------------------------------------|
| <b>en</b> Instructions for use     | <b>it</b> Naudojimo instrukcija       |
| <b>fr</b> Notice d'utilisation     | <b>lv</b> Lietošanas informācija      |
| <b>it</b> Istruzioni per l'uso     | <b>et</b> Kasutusjuhend               |
| <b>es</b> Instrucciones de uso     | <b>hu</b> Használati utasítás         |
| <b>pt</b> Instruções de utilização | <b>cs</b> Návod k použití             |
| <b>pl</b> Instrukcja użycia        | <b>sl</b> Navodila za uporabo         |
| <b>nl</b> Gebruiksaanwijzing       | <b>sk</b> Návod na používanie         |
| <b>sv</b> Bruksanvisning           | <b>bg</b> Инструкции за употреба      |
| <b>da</b> Brugsanvisningen         | <b>el</b> Οδηγίες χρήσης              |
| <b>no</b> Brugsanvisningen         | <b>tr</b> Kullanım kılavuzu           |
| <b>fi</b> Käyttöohje               | <b>uk</b> Інструкція для застосування |

# artPreOpaque<sup>PLUS</sup>

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# artPreOpaque PLUS

light-curing priming opaquer for use between  
artPrime PLUS and artOpaque PLUS.

## **en** Instructions for use, please read caefully!

### **Indication**

- Total or partial tooth loss
- Loss of hard tooth substance

### **Intended use**

Replacement of lost hard tooth substance by:

- Tooth-shaded fixture and for backing artVeneer®
- Veneering telescopic, attachment, bar and implant constructions with artVeneer®
- Tooth-shaded backing of trimmed acrylic teeth
- Fixing acrylic teeth to the model cast with a back protection plate or pontic
- Individualising acrylic teeth, artVeneer® and CAD/CAM restorations produced from M-PM® Disc and the artBloc® Temp
- Repairs

### **Contraindication**

- Hypersensitivity to one of the components

### **Components of the artConcept® PLUS system**

- artPrime PLUS - Metal bonding agent
- **artPreOpaque PLUS - Priming opaquer**
- artOpaque PLUS - Opaquer for covering framework
- artDentine PLUS/artEnamel PLUS Polymer - Tooth-shade PMMA acrylic
- artDentine Liquid NT PLUS (normal time) / artDentine Liquid LT PLUS (long time) - Monomer component for artDentine PLUS / artEnamel PLUS
- artConnect PLUS - Acrylic bonder
- artVeneer® - Tooth-shade PMMA anterior and posterior veneer

### **Composition**

Methacrylic ester

### **Processing**

- Metal primer (e.g. artPrime PLUS) is applied as a primer according to the instructions for use.
- To achieve optimal adhesion, artPreOpaque PLUS is applied using a brush with short, stiff bristles (no wash!), especially in the area of the retentions.
- Fluidity can be optimised by stirring on a mixing block.
- The material has very good light-curing properties that ensure optimal curing.

## Curing unit artPreOpaque PLUS

DC-XS / Uni-XS	Kulzer	90 sec
HeraFlash	Kulzer	90 sec
Labolight LV-II / III	GC	1 min
MPA 200	Dentsply	1 x HD
Solidilite EX / V	Shofu	1 min
Spektra 2000	Schütz Dental	3 min
Spektra 2000 fast	Schütz Dental	3 min
Spektra LED	Schütz Dental	30 sec
Spektramat	ivoclar vivadent	2 min
Triad II	Dentsply	5 min
Speed Labolight	Hager & Werken	1 min
Hilite power	Kulzer	90 sec

## Processing time of artPreOpaque PLUS

Depending on the ambient light 1- 3 min.

## Storage

Store in a dry and cool area at 10 - 25 °C / 50 - 77 °F without direct exposure to light.

## Disposal

The fully polymerized components of artConcept® PLUS are disposed of as residual waste. Disposal of non-polymerized product residues according to legal regulations. Do not allow entry into the sewage system or waterways.

## Notes

- Close all containers tightly after use.
- Do not use the material after the specified expiration date.
- To reduce the residual monomer content, the finished (sanded) but not polished polymerised workpiece should be stored in water at room temperature for about 48 hours before delivery.
- The residual monomer content after polymerisation of artDentine PLUS/artEnamel PLUS is < 1.8%.

## Danger and safety instructions (Contens <125ml)

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

---

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from [www.merz-dental.de](http://www.merz-dental.de).

---

**The product properties are based on compliance with and observation of these instructions for use.**

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**Merz Dental GmbH**