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[www.merz-dental.de](http://www.merz-dental.de)

ISO 13485

## (en) Instructions for use, please read carefully!

### Indications

Combitray® CC is distinguished by low shrinkage and high dimensional stability and is suitable for the fabrication of

- custom trays
- custom base plates
- registration sheets

### Composition

Combitray® CC powder contains

- PMMA
- PMMA-based terpolymer
- barbituric acid catalyst system
- inorganic pigments and fillers
- organic colorants

Combitray® CC liquid contains

- alkyl methacrylates
- alkyl dimethacrylate
- barbituric acid catalyst system
- silicone acrylate

### Dosage and application

- See table for mixing ratios and processing times.
- Remove the material and work it thoroughly by hand for a brief period.
- Roll out the material paste in sheets of even thickness using a rolling aid or similar, and then adapt.
- Cut the adapted material to size using a suitable instrument, e.g. a scalpel, and remove any excess material.
- The excess material can be used for modelling and for preparing a tray handle. For optimum bonding, the contact area can be moistened using a little Combitray® CC liquid.
- Following polymerisation, the tray is modified using a milling tool or similar.

### Times and quantities

Mixing ratio	9 ml : 25 g
Polymerisation time*	10 min*

- Mixing time\*
- Dwell time
- Processing time\* Adaption phase
- Max. processing time\* following adaption phase
- Polymerisation time\*

\*Based on a room temperature of 21 °C / 69,8 °F

## Trimming

To prevent inaccurate fitting after polymerisation, avoid excessive heat development during trimming and polishing.

## Special notes

- Close all containers tightly after use.
- Do not use the material after the specified expiration date.
- Store in a dry and cool area - up to 25 °C / 77 °F - away from direct sunlight.
- Record the LOT number with every process that requires identification of the material.

**The product properties are based on compliance with the instructions for use.**

---

## Combitray® CC powder

Contains PMMA, PMMA-based terpolymer, barbituric acid catalyst system, inorganic pigments and fillers, organic colorants

---

## Combitray® CC liquid

Contains MMA (methylmethacrylate), isobutyl methacrylate, silicone acrylate, ethylene glycol dimethacrylate

H225 Highly flammable liquid and vapour.  
H315 Causes skin irritation.  
H319 Causes serious eye irritation.  
H317 May cause an allergic skin reaction.  
H335 May cause respiratory irritation.  
P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.  
P280 Wear protective gloves/protective clothing/eye protection/face protection.  
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.  
P303+P361+P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.  
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.  
P337+P313 If eye irritation persists: Get medical advice/attention.  
P403+P233 Store in a well-ventilated place. Keep container tightly closed.  
P501 Disposal of contents/container in accordance with local/regional/national/international regulations.

Danger

---

For detailed information see the applicable Material Safety Data Sheet, which can be downloaded from [www.merz-dental.de](http://www.merz-dental.de)

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