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**en** Instructions for use, please read carefully!

Dental milling blanks made from zirconium dioxide

**Intended use**

For the manufacture of dental constructions incl. bridges ≥ 4 units with up to two connected pontics in the front and side tooth area.

**Indications**

• Loss of hard tooth substance

**Contraindications**

• Parafunctional habits  
• Inadequate preparation  
• Insufficient oral hygiene

• Insufficient space  
• Known intolerance to the contained components

**Processing/Notes**

Handling this medical device should only be allowed to trained staff. Shoulder or chamfer preparation is recommended. For tangential preparations, the crown border parameters in the construction design can be slightly increased. M-ZR white HT, M-ZR color HT and M-ZR multicolor ST is a sensitive high performance ceramic and should be processed with caution also in pre-sintered condition. Generally, the construction parameters need to be considered.

Please note: Depending on the construction, the connector cross section might need to be of larger dimensions. For example, in the case of wide-span bridges in the posterior region, the connector cross-section between two bridge pontics should be increased to at least 20 mm² if possible. An oval connector cross section is desirable; the height of the connector is the decisive factor. Frameworks for ceramic veneering need to be designed in a way to support the veneering ceramic in the cusp region and allow a constant layer thickness.

**Processing/Trimming**

Machining of blanks should always be conducted in the corresponding milling system. All information given by the manufacturer of the machine shall be taken into account. Information for Amann Girrbach users: Please find code (scale factor) for utilization of 98 mm discs in Amann Girrbach milling machines on blank margin. The frameworks need to be examined for visible defects. Do not process damaged or contaminated frameworks. The framework can be colored before final sintering M-ZR Basic Shade HC, M-ZR Pro Shade LC or M-ZR Art Shade (pay attention to separate instruction manual). The individual layer heights for design and positioning of the construction in the M-ZR multicolor ST blank can be taken from the following nesting table.

<b>Blank height (mm)</b>	<b>Layer 1+2: Incisal (mm / %)</b>	<b>Layer 3: Intermediate (mm / %)</b>	<b>Layer 4: Intermediate (mm / %)</b>	<b>Layer 5: Body / Dentine (mm / %)</b>
<b>14 mm</b>	3,5 / 24,9	2,1 / 15	2,1 / 15	6,3 / 45,1
<b>18 mm</b>	3,5 / 19,4	2,1 / 11,7	2,1 / 11,7	10,3 / 57,2
<b>22 mm</b>	3,5 / 15,9	2,1 / 9,6	2,1 / 9,6	14,3 / 64,9

**Please pay attention to our separate sintering instruction!**

**Sintering cycle for normal furnace filling without cover:**

↑ Heating up to 900 °C (1652 °F) with 8 °C/min (46 °F/min), → dwell at 900 °C (1652 °F) for 30 minutes, ↑ heating up to final temp. 1450 °C (2642 °F) with 3 °C/min (37 °F/min), → dwell at final sintering temperature 1450 °C (2642 °F) for 120 min, ↓ cooling to 200 °C (392 °F) with 10 °C/min (50 °F/min).

Avoid additional mechanical actions after sintering like blasting or grinding. If corrections are inevitable use water-cooled tools for the conditioning and make sure that there is no development of heat, which could lead to cracks in the material. Please work only with sharp, diamond-coated tools and very slight contact pressure. Areas that are under tension in clinical use (e.g. connectors) must not be processed. Do not separate in interdental regions. Avoid sharp edges.

**Processing / Ceramic veneering**

Please use a veneering ceramic with a suitable WAK and pay attention to the manufacturer's recommendation. Slowing down heating and cooling rate in case of massive constructions is strongly recommended.

<b>Weight per tooth unit [g]</b>	< 1	2	3	> 4
<b>Heating and cooling rate [°C/min]</b>	55	45	35	25
<b>Heating and cooling rate [°F/min]</b>	131	113	95	77

**Processing/Fixation**

For mounting we recommend conventional cementation with zinc oxide phosphates cements or glasionomer cements. Also luting composites may be used. Sufficient retentions and a stump with height of at least 3 mm are essential. A temporary mounting is not recommended!

**Storage**

The blanks are manufactured and tested to highest quality standards. Verify the integrity of the packaging and the blank before the first processing. Check if the content of the packaging corresponds to the declaration given on the label. Do not use damaged material. Storage only in the original containers in a cool and dry environment. Avoid vibrations, contaminations and contact with fluids.

**Disposal**

Residues from processing the M-ZR white HT, M-ZR color HT, M-ZR multicolor ST have to be disposed of as residual waste.

**Notes**

To ensure traceability and assignment of the product data as well as the LOT numbers for the patient at any time, it is recommended to store the milling blank in the original packaging during two milling operations.

**Warnings**

Processing zircon frameworks may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system on your milling machine and at the workplace is working properly for individual post-processing and that you are using personal protective equipment.

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from **www.merz-dental.de**.

**The product properties are based on compliance with and observation of these instructions for use.**

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