M-ZR multilayer HT

Zirconia Disc

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en Instructions for use

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en M-ZR multilayer HT – Instructions for Use

Safety notes

Please read the instructions for use carefully before taking the blank (made of zirconium dioxide) out of its packaging. They contain important information which is essential to ensure both the correct processing and the safety of patients and users. **M-ZR multilayer HT** blanks are manufactured and tested according to highest quality standards. In order to guaranty this level of quality during subsequent processing, the procedures described in the following must be absolutely adhered to.

General handling notes

Blanks made of **M-ZR multilayer HT** are supplied in a pre-sintered white-fired condition. They have limited stability and strength and a residual porosity. Therefore, these blanks must be handled with care. They must not be subjected to impact or vibration. Contamination must absolutely be avoided. Make sure restorations are only handled with dry and clean hands or gloves, and that they are under no circumstances contaminated with liquids (such as adhesives or marker pens). Coolants will reduce the material's translucency.

Indication:

M-ZR multilayer HT is a zirconium dioxide 5Y-TZP used for manufacturing fixed dental prostheses. The material is suitable for monolithic ceramic for prostheses involving partially or fully covered substructure for four or more units. For bruxism, the vertical dimension must be observed. In these cases, please consult with the dentist.

Contraindications:

In the event that there is an insufficient occlusal clearance and/or vertical prep wall, making the preparation unsuitable for an all-ceramic restoration, an alternative material must be chosen. Inlay bridges, endosseous implants and root posts are other contraindications.

Warning:

Known cross-reactions or interactions of this medical device with other medical products or material already present in the oral environment must be taken into consideration by the dentist when selecting this product.

Type/class of dental ceramic:

Type II/Class 5 (EN ISO 6872:2015)

Precautions:

Processing the material will result in the release of dust that may cause irritation of the eyes, skin or respiratory passages. Keep product dust away from eyes and avoid any contact with mucosa. After use, wash your hands. When using, do not eat or drink, do not ingest, and do not inhale dust particles. Use local suction and suitable mouth/ face protection during grinding.

Safety Data Sheets (SDS) may be obtained from <https://www.merz-dental.de> or by calling 0049 4381 403 400.

Adverse Effects:

We are not aware of any risks or adverse effects related to M-ZR multilayer HT Zirconia disc. If properly processed and used, adverse effects of this medical product are highly unlikely. However, when using zirconia, there is a potential for reactions of the immune system (such as allergies) or localized paresthesia (such as an irritating taste or irritation of the oral mucosa). If you are aware of an adverse effect or reaction, please contact Merz Dental and your dentist or physician immediately.

Operation instructions:

Caution:

The **M-ZR multilayer HT** is to be placed in direction of the arrowhead (upside = occlusal side .

Make sure that following criteria are secured for any restoration.

Thickness	
Anterior crown or bridge	≥ 0.8 mm
Posterior crown or bridge, inlay or onlay	≥ 1.0 mm
Veneer	≥ 0.4 mm
Connector cross-section	
Anterior 2- or 3-unit bridges	≥ 12 mm ²
Posterior 2- or 3-unit bridges	≥ 16 mm ²
Unit	
Bridges involving molar restoration	≤ 3 units
Bridges NOT involving molar restoration	≤ 6 units

Do not use this product for fabrication of bridge restorations including 2- or more connecting pontics.

When building up porcelain on a frame, build up and shape individually.

Fix a disc according to operating instructions of the CAD/CAM system and enter the set value (see .

Start the system to cut and machine the disc.

Neither coolant nor compressed air should be used during processing. Please observe the instructions for use of your milling machine as well as the parameters of the CAD/CAM software.

Please separate the finished milled parts from the blanks carefully using suitable tools. Then, thin the edges which have been thickened by the machine and remove the tapping.

Visual inspection:


Before the milled mounts are processed further, they must be inspected for the following faults: - Shiny areas on the surface (indicating a worm milling cutter) - Discoloration - Material spalling (due to the milling strategy and milling cutter) - Cracks - Faulty restorations must not be processed further.

Final Sintering:

Sintering can take place in all common dental sintering furnaces which are approved for the sintering of restorations made of zirconium dioxide. Since zirconium dioxide is known to be a poor thermal conductor, it is recommended that you slowly heat the mounts to the required temperature (see sintering graph) and cool them down equally slowly.

Recommended Sinter-Program:

Increase temperature by 5 °C/minute from room temperature to 1450 °C. Keep at 1450 °C for 120 minutes, and then cool from 1450 °C to room temperature decreasing the temperature by 5 °C/minute.

Further processing of the finally-sintered restorations:

IMPORTANT: The sintered units will only achieve their final color after glaze firing. If possible, the -sintered units should not be milled by hand. However, if this is unavoidable, only water-cooled diamond-fitted tools in correct working order may be used. Ideally, the interdental connection points should not be polished. In principle, the basal polishing of these points must always be avoided for stability reasons (predetermined breaking points). For abutments, sharp edges should be avoided and rounded, if possible.

Grinding:

In order to protect the antagonists (abrasion) and for material engineering reasons, ground occlusal contact points and surfaces must be either high gloss polished after the try-in and/or glaze-fired.

Cementation:

The interior surfaces of the restoration should be receiving mechanical retention by means of blasting with 29.0-43.5 psi (0.2-0.3 Mpa) -Particle size of Al₂O₃, 50-100 µm. Before the cement is applied internally, the blasted interior surface should be cleaned with alcohol. State-of-the-art dental technology prefers self-adhesive and adhesive cementation.

Material:	ZrO ₂ Y-TZP
Dimensions:	98 mm Ø X 14/18/22
Shades:	<ul style="list-style-type: none"> • M-ZR multilayer HT – A1 (5 - Layer multilayer disc) • M-ZR multilayer HT – A2 (5 - Layer multilayer disc) • M-ZR multilayer HT – A3 (5 - Layer multilayer disc)

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MT Promedt Consulting GmbH Altenhofstrasse 80 66386 St. Ingbert, Germany