PremEco® Line Gießkunststoff

Kaltpolymerisat

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PremEco® Line casting resin

Cold-curing polymer

USA Instructions for use, please read carefully!

Caution: Federal law restricts this device to sale by or on the order of a licensed dental practitioner.

Indications

Colour-fast, cold-curing, customisable methylmethacrylate acrylic base specially developed for casting partial and full dentures.

Classification according to MDD 93/42/EEC and MDR (EU) 2017/745

Class IIa medical device for long-term contact with intact mucous membranes.

Patient target groups

All genders and age groups

Intended users

Dental professionals

Composition

Polymer	Monomer
PMMA (polymer and copolymer)	MMA (methylmethacrylate)
barbituric acid catalyst system	dimethacrylate
organic colorants	barbituric acid catalyst system
inorganic pigments	

Dosage and application


- See table for mixing ratios and processing times.
- Use class IV plaster for making casts where possible.
- Always wet casts sufficiently for 5 to 15 minutes depending on the degree of dryness.
- Use standard alginate isolation to separate plaster from acrylic.
- Do not compress the powder when using dosage aids.

Times and quantities

Mixing ratio	30 g : 18 ml
Polymerisation temperature	45 °C / 113 °F
Polymerisation time	15 min ¹ / 30 min ²
Pressure required during polymerisation	2 - 2,5 bar

¹ in matrix technique, plaster, putty

² in flask duplicating technique, duplicating gel, duplicating silicone

PremEco® Line casting resin is a product of the PremEco® Line production system. For more information see www.merz-dental.de.

Trimming

To prevent inaccurate fitting after polymerisation, avoid excessive heat development during trimming and polishing.

Bond

- Roughen the basal surfaces of the acrylic teeth before bonding.
- Ensure that no impurities such as wax or isolating agent are present.
- Mechanical retention is recommended.
- Adequate mechanical retention is essential for ceramic teeth.

Restoration

Products of PremEco® Line casting resin can be restored and supplemented at any time with all standard MMA-based cold-curing polymers, such as Weropress®.

Special notes

- Close all powder and liquid containers tightly after use.
- Store in a dry and cool area away from direct sunlight.
- Take precautionary measures against static discharges.
- Do not use the material after the expiration date.
- Record the LOT number with every process that requires identification of the material.
- The residual monomer content after polymerisation of PremEco® Line casting resin is < 2,5%.
- For more information on processing acrylics and reducing residual monomer content, please visit www.merz-dental.de/service
- All serious incidents occurring with users and / or patients in connection with this product must be reported to Merz Dental, as well as to the German Federal Institute for Drugs and Medical Devices (BfArM), or the respective competent authority of the member state in which the user and / or patient is resident.

Disposal information

Disposal of contents/container in accordance with local/regional/national/international regulations.

Danger and safety instruction PremEco® Line casting resin Monomer

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P233 Keep container tightly closed.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P262 Do not get in eyes, on skin, or on clothing.

P210 Keep away from heat. No Smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

Danger and safety instruction PremEco® Line casting resin Polymer

May cause an allergic skin reaction.

Avoid breathing dust/fume/gas/mist/vapours/spray.

Classification

PremEco® Line casting resin conforms to EN ISO 20795-1 type 2 class 2.

The product properties are based on compliance with the instructions for use.

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de.

Date of information 2022-12

PremEco® Line casting resin

Cold-curing polymer

(en) Instructions for use, please read carefully!

Intended purpose

Colour-fast, cold-curing, customisable methylmethacrylate acrylic base specially developed for casting partial and full dentures.

Classification according to MDD 93/42/EEC and MDR (EU) 2017/745

Class IIa medical device for long-term contact with intact mucous membranes.

Patient target groups

All genders and age groups

Intended users

Dental professionals

Indication

- Total or partial tooth loss

Contraindication

- Hypersensitivity to one of the components

Composition

Polymer	Monomer
PMMA (polymer and copolymer)	MMA (methylmethacrylate)
barbituric acid catalyst system	dimethacrylate
organic colorants	barbituric acid catalyst system
inorganic pigments	

Dosage and application


- See table for mixing ratios and processing times.
- Use class IV plaster for making casts where possible.
- Always wet casts sufficiently for 5 to 15 minutes depending on the degree of dryness.
- Use standard alginate isolation to separate plaster from acrylic.
- Do not compress the powder when using dosage aids.

Times and quantities

Mixing ratio	30 g : 18 ml
Polymerisation temperature	45 °C / 113 °F
Polymerisation time	15 min ¹ / 30 min ²
Pressure required during polymerisation	2 - 2,5 bar

¹ in matrix technique, plaster, putty

² in flask duplicating technique, duplicating gel, duplicating silicone

PremEco® Line casting resin is a product of the PremEco® Line production system. For more information see www.merz-dental.de.

Trimming

To prevent inaccurate fitting after polymerisation, avoid excessive heat development during trimming and polishing.

Bond

- Roughen the basal surfaces of the acrylic teeth before bonding.
- Ensure that no impurities such as wax or isolating agent are present.
- Mechanical retention is recommended.
- Adequate mechanical retention is essential for ceramic teeth.
- Use artConnect PLUS bonding agent to support bonding.

Restoration

Products of PremEco® Line casting resin can also be restored and supplemented at any time with all standard MMA-based cold-curing polymers. Following roughening of the surface (sanding / blasting / grinding), further conditioning is recommended, preferably with artConnect PLUS or a similar MMA-based bonding agent.

Storage

- Store in a dry and cool area
- Avoid direct sunlight

Special notes

- Close all powder and liquid containers tightly after use.
- Changes to the construction through use, e.g. due to fissures, cracks, clouding or significant signs of abrasion, necessitate a repair or refabrication. Monitoring should be carried out with regular assessment of the denture at the follow-up presentation in the dental surgery.
- Do not use the material after the expiration date.
- Record the LOT number with every process that requires identification of the material.
- The residual monomer content after polymerisation of PremEco® Line casting resin is < 2,5%.
- Processing frameworks may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system on your milling machine and at the workplace is working properly for individual post-processing and that you are using personal protective equipment.
- For more information on processing acrylics and reducing residual monomer content, please visit www.merz-dental.de/service
- All serious incidents occurring with users and / or patients in connection with this product must be reported to Merz Dental, as well as to the German Federal Institute for Drugs and Medical Devices (BfArM), or the respective competent authority of the member state in which the user and / or patient is resident.

Disposal

- The fully polymerized components of PremEco Line casting resins are disposed of as residual waste.
- Disposal of nonpolymerized product residues according to legal regulations.
- Do not allow entry into the sewage system or waterways.

Danger and safety instruction PremEco® Line casting resin Monomer

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P233 Keep container tightly closed.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P262 Do not get in eyes, on skin, or on clothing.

P210 Keep away from heat. No Smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

Danger and safety instruction PremEco® Line casting resin Polymer

May cause an allergic skin reaction.

Avoid breathing dust/fume/gas/mist/vapours/spray.

Classification

PremEco® Line casting resin conforms to EN ISO 20795-1 type 2 class 2.

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de.

The product properties are based on compliance with the instructions for use.

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Merz Dental GmbH