artDentine**PLUS**

artEnamel**PLUS**

artDentine**PLUS** Liquid NT

artDentine**PLUS** Liquid LT

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- |                               |                             |
|-------------------------------|-----------------------------|
| (en) Instructions for use     | (it) Naudojimo instrukcija  |
| (fr) Notice d'utilisation     | (lv) Lietošanas informācija |
| (it) Istruzioni per l'uso     | (et) Kasutusjuhend          |
| (es) Instrucciones de uso     | (hu) Használati utasítás    |
| (pt) Instruções de utilização | (cs) Návod k použití        |
| (pl) Instrukcja użycia        | (sl) Navodila za uporabo    |
| (nl) Gebruiksaanwijzing       | (sk) Návod na používanie    |
| (sv) Bruksanvisning           | (bg) Инструкции за употреба |
| (da) Brugsanvisningen         | (el) Οδηγίες χρήσης         |
| (no) Brugsanvisningen         | (tr) Kullanım kılavuzu      |
| (fi) Käyttöohje               |                             |

# artDentine/Enamel**PLUS** artDentine**PLUS** Liquid NT/LT

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Hersteller

# artDentine/Enamel<sup>PLUS</sup> artDentine<sup>PLUS</sup> Liquid NT/LT Polymer and Monomer

Tooth-shaded auto-polymerizing PMMA acrylic, shade-matched to artVeneer® / artVeneer® life and Merz Dental acrylic teeth.

## **(en) Instructions for use, please read caefully!**

### **Indication**

- Total or partial tooth loss
- Loss of hard tooth substance

### **Intended purpose**

Replacement of lost hard tooth substance by:

- Tooth-shaded fixture and for backing artVeneer® / artVeneer® life
- Veneering telescopic, attachment, bar and implant constructions with artVeneer® / artVeneer® life
- Tooth-shaded backing of trimmed acrylic teeth
- Fixing acrylic teeth to the model cast with a back protection plate or pontic
- Individualising acrylic teeth, artVeneer® / artVeneer® life and CAD/CAM restorations produced from M-PM® Disc and the artBloc® Temp
- Repairs

### **Classification according to MDD 93/42/EEC and MDR (EU) 2017/745**

Class IIa medical device for long-term contact with intact mucous membranes.

### **Patient target groups**

All genders and age groups

### **Intended users**

Dental professionals

### **Contraindication**

- Hypersensitivity to one of the components

### **Components of the artConcept® PLUS system**

- artPrime PLUS - Metal bonding agent
- artPreOpaque PLUS - Priming opaquer
- artOpaque PLUS - Opaquer for covering framework
- artDentine PLUS/artEnamel PLUS Polymer - Tooth-shade PMMA acrylic
- artDentine PLUS Liquid NT (normal time) / artDentine PLUS Liquid LT (long time) - Monomer component for artDentine PLUS / artEnamel PLUS
- artConnect PLUS - Acrylic bonder
- artVeneer® - Tooth-shade PMMA anterior and posterior veneer

### **Composition**

artDentine / artEnamel PLUS	artDentine PLUS Liquid NT	artDentine PLUS Liquid LT
PMMA copolymer	MMA (methylmethacrylate)	MMA (methylmethacrylate)
barbituric acid derivative	dimethacrylate	ethermethacrylate
organic colourants	barbituric acid catalyst system	barbituric acid catalyst system
inorganic colourants		

## Processing

### Fixture of the veneers

- Remove remnants of wax and fat, grinding dust and other contamination from the artVeneer® / artVeneer® life veneers.
- After evaporation, ensure complete drying, also in the undercuts. There must be no contamination from blowing off with compressed air (e.g. from oil residues in the compressed air system).
- When roughening artVeneer® / artVeneer® life or acrylic teeth, affix small groove retentions as required for larger components.
- Wetten the surface to be joined with artConnect PLUS and apply sufficient that the connecting surface remains moist for 5 min to work in. If the connecting surface with artConnect PLUS dries out during this time, moisten again with artConnect PLUS. The exposure time must not be less than 5 min; the fluid artDentine PLUS / artEnamel PLUS resin has to be applied to the conditioned surface after 10 min at the latest.
- artEnamel PLUS can be used for individualisation in the incisal area.

**Tip:** To increase the translucency, artDentine PLUS polymer can be mixed individually with artEnamel PLUS (TC Transpa Clear or artEnamel PLUS incisal).

### Dosage

	artDentine PLUS/artEnamel PLUS with artDentine PLUS Liquid NT (normal time)	artDentine PLUS/artEnamel PLUS with artDentine PLUS Liquid LT (long time)
Mixing ratio	Polymer: Liquid, 10 g : 7 ml	Polymer: Liquid, 10 g : 7g
Mixing time	30 sec	30 sec
Dwelling time	65 sec*	155 sec* / 50 sec**
Fluid phase	90 sec*	190 sec* / 100 sec**
Polymerisation	15 min in the pressure pot at 45 °C / 113 °F and 2 - 2,5 bar pressure	in the pressure pot at 45 °C / 113 °F and 2 - 2,5 bar pressure

#### artDentine PLUS/artEnamel PLUS with artDentine PLUS Liquid NT (normal time)

<b>30 sec</b>	<b>65 sec*</b>	<b>90 sec*</b>	<b>15 min</b>
Mixing time	Dwelling time	Fluid phase	Polymerisation

#### artDentine PLUS/artEnamel PLUS with artDentine PLUS Liquid LT (long time)

<b>30 sec</b>	<b>155 sec* / 50 sec**</b>	<b>190 sec* / 100 sec**</b>	<b>20 min</b>
Mixing time	Dwelling time	Fluid phase	Polymerisation

\*Based on a room temperature of 21 °C / 69,8 °F

\*\*Based on a room temperature of 30 °C / 86 °F

Use class IV plaster for making casts where possible. Always sufficiently wetten the models and use standard alginate isolation to separate plaster from acrylic.

### Trimming

- Cross-cut carbide burs suitable for resin or suitable cutting discs should be used to rework.
- To prevent inaccurate fitting after polymerisation, avoid excessive heat development during trimming and polishing.

### Restoration

Products of artDentine PLUS can also be restored and supplemented at any time with all standard MMA-based cold-curing polymers. Prior conditioning of the surface is necessary, preferably with artConnect PLUS or a similar acrylate-based bonding agent in accordance with the specifications of the manufacturer.

### Storage

Store in a dry and cool area at 10 - 25 °C / 50 - 77 °F without direct exposure to light.

## Disposal

The fully polymerized components of artConcept® PLUS are disposed of as residual waste. Disposal of non-polymerized product residues according to legal regulations. Do not allow entry into the sewage system or waterways.

## Notes

- Close all containers tightly after use.
- Do not use the material after the specified expiration date.
- All serious incidents occurring with users and / or patients in connection with this product must be reported to Merz Dental, as well as to the German Federal Institute for Drugs and Medical Devices (BfArM), or the respective competent authority of the member state in which the user and / or patient is resident.
- The residual monomer content after polymerisation of artDentine PLUS/artEnamel PLUS is < 1.8%
- For more information on processing acrylics and reducing residual monomer content, please visit [www.merz-dental.de/service](http://www.merz-dental.de/service).

## Danger and safety instructions

### artDentine PLUS Liquid NT/LT

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P210 Keep away from heat. No Smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P403+P233 Store in a well-ventilated place. Keep container tightly closed.

---

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from [www.merz-dental.de](http://www.merz-dental.de).

---

**The product properties are based on compliance with and observation of these instructions for use.**

Date of information 2021-07

## Symbols

 Article number

Protect from the sun

 Batch

Observe Instructions for use (IfU)

 Use by

Temperature limitation

 Reference that the product is a medical device

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**Merz Dental GmbH**