artOpaque^{PLUS} GUM

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- | | |
|------------------------------------|---------------------------------------|
| en Instructions for use | lt Naudojimo instrukcija |
| fr Notice d'utilisation | lv Lietošanas informācija |
| it Istruzioni per l'uso | et Kasutusjuhend |
| es Instrucciones de uso | hu Használati utasítás |
| pt Instruções de utilização | cs Návod k použití |
| pl Instrukcja użycia | sl Navodila za uporabo |
| nl Gebruiksaanwijzing | sk Návod na používanie |
| sv Bruksanvisning | bg Инструкции за употреба |
| da Brugsanvisningen | el Οδηγίες χρήσης |
| no Bruksanvisningen | tr Kullanım kılavuzu |
| fi Käyttöohje | uk Інструкція для застосування |

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(en) Instructions for use, please read carefully!

Intended purpose

Framework covering materials are composite-based materials for covering different framework materials in preparation for subsequent veneering with veneering composite or denture resin.

General description of product

- Light-cured single-component material for colour coating a framework material made from standard commercial dental metals and alloys. The material is light-cured and is polymerised in the light-curing unit. artOpaque PLUS GUM also contains an adhesive component that guarantees bonding to the framework material.
Thus it is not necessary to use an additional adhesive bonding system.
- Before completion with plastic, the material is applied to the framework retentions to form a gingiva-coloured coating.
- artOpaque PLUS GUM is available in GUM light and GUM dark shades.
Both shades can be mixed to form customised hues.

Users

For use by laboratory technicians in a dental laboratory.

Composition

- Silicon dioxide, glass powder, 2-[[[(butylamino)carbonyl] oxy] ethyl acrylate, di-urethane dimethacrylate, pigments, initiators
- Fillers: 45 wt.%, inorganic filling materials (0.005 - 3.0 µm)

Indications

For coating model cast retentions (standard commercial dental metals and alloys) and gingival finishing.

Contraindications

Do not use in the case of a known allergy to one of the components.

Warnings

Contains: 2-[[[(butylamino)carbonyl] oxy] ethyl acrylate, di-urethane dimethacrylate diphenyl (2,4,6- trimethylbenzoyl) phosphine oxide. May cause an allergic skin reaction.

Safety instructions

Avoid breathing vapours/ spray. Wear protective gloves.

Interactions with other medicinal products

Phenolic substances such as Eugenol inhibit polymerisation. Therefore, do not use any material containing these substances.

Application / Preparation

Processing times: 1 - 3 minutes, depending on lighting conditions.

A. Preparatory work

The model casting framework is produced in accordance with the manufacturer's specifications. The surface is sand-blasted with aluminium oxide 50-110 µm at 2 bar and blown off with oil-free compressed air.

B. Procedure

A thin (wash-type) layer of the artOpaque *PLUS* GUM is applied with a stiff short-hair brush. At least two layers must be applied to ensure that the coating covers the whole surface. A transparent opaquer layer has a negative effect on the colour.

The flowability of the opaquer can be optimised by stirring it on a mixing pad. Ensure that the applied layer is extra thin around the retentions. Each layer is polymerised separately. (see Polymerisation Times table)

Light-curing unit	Time
Spektra LED	1 min
Spektra 2000	3 min
HiLite / UniXS	90 sec
Labolight LV-II / III	1 min
Solidilite	1 min

C. Subsequent processing

The plastic can then be cast on.

Notes on processing

A light polymerisation unit with an emission spectrum of at least 310-500 nm should be used. The physical characteristics required can only be achieved if the correct lamps are used. Regular checks of the light intensity in accordance with the manufacturer's specifications are therefore required. Do not use intermediate polymerisation units.

Troubleshooting / FAQ List

- Remove any layers of grease and polishing waste from the frame surfaces.
- Do not mix artOpaque *PLUS* GUM with other opaquer liquids / powders.
- Opaquer layers that have been applied too thickly prevent optimal polymerisation and thus weaken the bond. Curing cannot be improved by extending the polymerisation time.

Fault	Cause	Corrective action
does not solidify	applied too thickly	only applied in a wash-like coating on multiple layers
	inadequate polymerisation	observe the polymerisation times - check lamp and replace if necessary - curing cannot be improved by extending the polymerisation time - do not use intermediate polymerisation lamps
opaquer peels off	opaquer not properly polymerised / applied too thickly	only applied in a wash-like coating or in very thin layers
		wrong polymerisation times / check lamp and replace if necessary
surface greasy	inadequate polymerisation	observe the polymerisation times
		check device / service device regularly

Storage

- storage temperature 10-25°C / 50-77°F
- Close the syringe tightly immediately after use.

Shelf life

The maximum shelf life is printed on the label of each syringe. Do not use after the expiry date.

Side effects

With proper preparation and use of this medical device, adverse effects are extremely rare. However, immune reactions (such as allergies) or local discomfort cannot in principle be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case.

Disposal

Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

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MANI MEDICAL GERMANY GmbH
Merz Dental GmbH