# artOpaque<sup>PLUS</sup>

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- |                               |                                  |
|-------------------------------|----------------------------------|
| (en) Instructions for use     | (lt) Naudojimo instrukcija       |
| (fr) Notice d'utilisation     | (lv) Lietošanas informācija      |
| (it) Istruzioni per l'uso     | (et) Kasutusjuhend               |
| (es) Instrucciones de uso     | (hu) Használati utasítás         |
| (pt) Instruções de utilização | (cs) Návod k použití             |
| (pl) Instrukcja użycia        | (sl) Navodila za uporabo         |
| (nl) Gebruiksaanwijzing       | (sk) Návod na používanie         |
| (sv) Bruksanvisning           | (bg) Инструкции за употреба      |
| (da) Brugsanvisningen         | (el) Οδηγίες χρήσης              |
| (no) Bruksanvisningen         | (tr) Kullanım kılavuzu           |
| (fi) Käyttöohje               | (uk) Інструкція для застосування |

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Stand der Informationen 2024-02

**(en) Instructions for use, please read carefully!**

## **Intended purpose**

Framework covering materials are composite-based materials for covering different framework materials in preparation for subsequent veneering with veneering composite or denture resin

## **General description of product**

Light-curing paste opaquers for crown and bridge techniques that are suitable for use with composite veneers.

## **Users**

For use by laboratory technicians in a dental laboratory.

## **Composition**

- Silicon dioxide, glass powder, 2-[[butylamino)carbonyl] oxy] ethyl acrylate, di-urethane dimethacrylate, pigments, initiators
- Fillers: 49 wt.%, inorganic filling materials (0.002-3.0 µm)

## **Indications**

For coating framework materials used for crown and bridge techniques (standard commercial dental metals and alloys) for new items and repairs.

## **Contraindications**

If the patient is allergic or hypersensitive to one of the components, this product must not be used or only under the strict supervision of the treating doctor/dentist.

## **Warnings**

Contains: 2-[[butylamino)carbonyl] oxy] ethyl acrylate, di-urethane dimethacrylate diphenyl (2,4,6- trimethylbenzoyl) phosphine oxide.

May cause an allergic skin reaction.

## **Safety instructions**

Avoid breathing vapours/ spray. Wear protective gloves.

## **Interactions with other medicinal products**

Phenolic substances such as Eugenol inhibit polymerisation. Therefore, do not use any material containing these substances.

## **Application / Preparation**

Processing times: 1-3 minutes, depending on lighting conditions.

## A. Preparatory work

- The frameworks are modelled, cast or CAD/CAM milled and finished as usual in accordance with the applicable guidelines on dental procedures.
- Standard commercial dental metals and alloys, PMMA, PEEK and zirconia may be used as framework materials.
- All materials must be processed in accordance with the manufacturer's specifications and prepared for further processing (sand-blasting, conditioning).
- The frameworks must be conditioned before the opaquer is used. This conditioning varies depending on the framework material used. Observe the information provided by the manufacturer.

## B. Procedure

A thin (wash-type) layer of the artOpaque *PLUS* is applied with a stiff short-hair brush. At least two layers must be applied to ensure that the coating covers the whole surface. A transparent opaquer layer has a negative effect on the colour.

The flowability of the opaquer can be optimised by stirring it on a mixing pad. Ensure that the applied layer is extra thin around the retentions. Each layer is polymerised separately. (see Polymerisation Times table)

| Light-curing unit     | Time   |
|-----------------------|--------|
| Spektra LED           | 1 min  |
| Spektra 2000          | 3 min  |
| HiLite / UniXS        | 90 sec |
| Labolight LV-II / III | 1 min  |
| Solidilite            | 1 min  |

## C. Subsequent processing

The restoration is then designed with a commercially available veneering composite and PMMA-based materials. Observe the information provided by the manufacturer.

## D. Repairs

- All repairs are extraoral.
- The surface to be repaired is prepared as described in A.
- The procedure is described in B and subsequent processing in C.

## Troubleshooting / FAQ List

- Remove any layers of grease and polishing waste from the frame surfaces.
- Do not mix artOpaques *PLUS* with other opaquer liquids / powders.
- Opaquer layers that have been applied too thickly prevent optimal polymerisation and thus weaken the bond. Curing cannot be improved by extending the polymerisation time.
- Do not use intermediate polymerisation units for polymerisation.

| <b>Fault</b>      | <b>Cause</b>   | <b>Corrective action</b>  |
|-------------------|--|---|
| does not solidify | applied too thickly                                    | only applied in a wash-like coating on multiple layers  |
|                   | inadequate polymerisation                              | observe the polymerisation times<br>- check lamp and replace if necessary<br>- curing cannot be improved by extending the polymerisation time<br>- do not use intermediate polymerisation lamps |
| opaquer peels off | opaquer not properly polymerised / applied too thickly | only applied in a wash-like coating or in very thin layers  |
|                   |  | wrong polymerisation times / check lamp and replace if necessary  |
| surface greasy    | inadequate polymerisation                              | observe the polymerisation times  |
|                   |  | check device / service device regularly   |

### **Storage**

- storage temperature 10-25°C / 50-77°F
- Close syringe carefully.

### **Shelf life**

The maximum shelf life is printed on the label of each syringe. Do not use after the expiry date.

### **Side effects**

With proper preparation and use of this medical device, adverse effects are extremely rare. However, immune reactions (such as allergies) or local discomfort cannot in principle be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case.

### **Disposal**

Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

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**MANI MEDICAL GERMANY GmbH**
**Merz Dental GmbH**