# artPreOpaque<sup>PLUS</sup>

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- |                                    |                                       |
|------------------------------------|---------------------------------------|
| <b>en</b> Instructions for use     | <b>lt</b> Naudojimo instrukcija       |
| <b>fr</b> Notice d'utilisation     | <b>lv</b> Lietošanas informācija      |
| <b>it</b> Istruzioni per l'uso     | <b>et</b> Kasutusjuhend               |
| <b>es</b> Instrucciones de uso     | <b>hu</b> Használati utasítás         |
| <b>pt</b> Instruções de utilização | <b>cs</b> Návod k použití             |
| <b>pl</b> Instrukcja użycia        | <b>sl</b> Navodila za uporabo         |
| <b>nl</b> Gebruiksaanwijzing       | <b>sk</b> Návod na používanie         |
| <b>sv</b> Bruksanvisning           | <b>bg</b> Инструкции за употреба      |
| <b>da</b> Brugsanvisningen         | <b>el</b> Οδηγίες χρήσης              |
| <b>no</b> Bruksanvisningen         | <b>tr</b> Kullanım kılavuzu           |
| <b>fi</b> Käyttöohje               | <b>uk</b> Інструкція для застосування |

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**(en) Instructions for use, please read carefully!**

## **Intended purpose**

Framework covering materials are composite-based materials for covering different framework materials in preparation for subsequent veneering with veneering composite or denture resin

## **General description of product**

Light-cured universal priming opaquer for use between metal primer and major commercially available light-cured powder/paste opaquers.  
artPreOpaque PLUS improves the chemical bond between primer and opaquer and contributes to easy, safe and long-lasting adhesive bonding.

## **Users**

For use by laboratory technicians in a dental laboratory.

## **Composition**

- Silicon dioxide, di-urethane dimethacrylate, 2-[[[(butylamino)carbonyl]oxy]ethyl acrylate, glass powder, pigments, initiators
- Fillers: 62 wt.%, inorganic filling materials (0.005-3.0 µm)

## **Indications**

For coating framework materials used for crown and bridge techniques (standard commercial dental metals and alloys) for new items and repairs.

## **Contraindications**

Do not use in the case of a known allergy to one of the components.

## **Warnings**

Warning. Contains di-urethane dimethacrylate, 2-[[[(butylamino)carbonyl]oxy]ethyl acrylate, diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide.  
May cause allergic skin reaction.

## **Safety instructions**

Avoid breathing vapours/ spray. Wear protective gloves.

## **Interactions with other medicinal products**

- Phenolic substances such as Eugenol inhibit polymerisation. Therefore, do not use any material containing these substances.
- The dentist should consider known interactions and cross-reactions of the medical device with other materials already in the patient's mouth before using the product.

## A. Preparatory work

- The frameworks are modelled, cast or CAD/CAM milled and finished as usual in accordance with the applicable guidelines on dental procedures.
- Standard commercial dental metals and alloys, PMMA, PEEK and zirconia may be used as framework materials.
- All materials must be processed in accordance with the manufacturer's specifications and prepared for further processing (sand-blasting, conditioning).
- The frameworks must be conditioned before the opaquer is used. This conditioning varies depending on the framework material used. Observe the information provided by the manufacturer.

## B. Procedure

A thin layer of artPreOpaque *PLUS* is applied with a stiff short-hair brush. At least two layers must be applied to ensure optimal adhesive bonding.

The flowability of the opaquer can be optimised by stirring it on a mixing pad. Ensure that the applied layer is extra thin around the retentions. Each layer is polymerised separately. (see Polymerisation Times table)

| Light-curing unit     | Time   |
|-----------------------|--------|
| HiLite / UniXS        | 90 sec |
| Labolight LV-II / III | 1 min  |
| Solidilite            | 1 min  |

## C. Subsequent processing

A standard commercial opaquer is then applied and a standard commercial composite veneer used for forming. Observe the information provided by the manufacturer.

## D. Repairs

- All repairs are extraoral.
- The surface to be repaired is prepared as described in A.
- The procedure is described in B and subsequent processing in C.

## Troubleshooting / FAQ List

- Remove any layers of grease and polishing waste from the frame surfaces.
- Do not use intermediate polymerisation units for polymerisation.

| <b>Fault</b>                         | <b>Cause</b>               | <b>Corrective action</b>  |
|--------------------------------------|----------------------------|---|
| does not solidify                    | applied too thickly        | only applied on multiple layers   |
|                                      | inadequate polymerisation  | observe the polymerisation times<br>- check lamp and replace if necessary<br>- curing cannot be improved by extending the polymerisation time<br>- do not use intermediate polymerisation lamps |
| material is too solid in the syringe | syringe stored below 10 °C | store material at room temperature or bring to room temperature in good time prior to processing  |
| surface greasy                       | inadequate polymerisation  | observe the polymerisation times  |
|                                      |                            | check device / service device regularly   |

### **Fabrication of frameworks and preparation**

Processing times: 1-3 minutes, depending on lighting conditions.

### **Storage**

- Storage temperature 10-25°C / 50-77°F
- Close the syringe tightly immediately after use.

### **Shelf life**

The maximum shelf life is printed on the label of each syringe. Do not use after the expiry date.

### **Side effects**

With proper preparation and use of this medical device, adverse effects are extremely rare. However, immune reactions (such as allergies) or local discomfort cannot in principle be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case.

### **Disposal**

Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

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**MANI MEDICAL GERMANY GmbH**
**Merz Dental GmbH**