# PremEco® Line

# Präzisions

# Gel

## Gebrauchsinformation

- (en) Instructions for use
- (fr) Mode d'emploi
- (it) Istruzioni per l'uso
- (es) Instrucciones de uso
- (pt) Instruções de utilização
- (pl) Instrukcja użycia
- (nl) Gebruiksaanwijzing
- (sv) Information för användning
- (da) Brugsanvisningen
- (no) Brugsanvisningen
- (fi) Käyttöohje
- (lt) Informacinis lapelis
- (lv) Lietošanas informācija
- (et) Kasutusteave
- (hu) Használati tájékoztató
- (cs) Informace k použití
- (sl) Navodila za uporabo
- (sk) Návod na použitie
- (el) Οδηγίες χρήσης
- (tr) Kullanım bilgidir

DENTAL

# PremEco® Line Präzisionsgel

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# PremEco® Line precision gel

Agar-agar base

**en** Instructions for use, please read carefully!

## Indications

Particularly suitable for manufacturing plastic dentures using pouring method and for production of embedded models for model casting.

## Composition

PremEco® Line precision gel contains

- water, glycerine, hydrocolloid
- colorant, odorant and preservatives

## Dosage and application

- The gel supplied in the bucket is ready to be used.
- Cut the solid gel into smaller pieces before putting it into the duplicating unit to allow it to melt faster.
- Melting temperature 95 °C
- Processing temperature between 47 °C and 49 °C
- The gel can be used up to 10 times depending on the contamination.
- If the gel has been used previously, remove residues and impurities (embedding compound, plaster etc.) before melting it again.
- We recommend adding approx. 5% water every time the gel is melted, because water can evaporate during the process.
- Do not overfill the duplicating unit (no more than 80%) to prevent overflow when boiling.

**TIP:** ensure that the ventilation operates correctly.

- Check the temperature of the gel during the melting process (95 °C) and during processing (47-49 °C) occasionally with an external thermometer, because the temperature must be maintained exactly for the best results. If any deviations in the temperature setting on the unit are found, the unit must be calibrated.
- Pour the gel into the flask.
- Leave the duplicate to cool for 45 minutes in a 3-5 cm high water bath.

## Melting in the microwave

- Melt approximately 350 g gel cut into small pieces with 20 ml water in a suitable vessel in the microwave for approximately three minutes, depending on the power of the unit.
- The gel must not boil, if necessary interrupt the melting process and stir the gel briefly.
- After melting the gel, allow it to cool to 47 °C to 49 °C while stirring.

PremEco® Line precision gel is a product of the PremEco® Line production system.

For more information see [www.merz-dental.de](http://www.merz-dental.de)

### Special notes

- Always store new and used material separately to prevent contamination of new material.
- The gel can no longer be used after the tenth melting process and it loses its properties. It can be disposed of in the household garbage.
- Gel duplicating material is not suitable for making plaster duplicates.
- Close PremEco® Line precision gel containers tightly after use to prevent the gel from drying out.
- Store in a dry and cool area away from direct sunlight.
- Do not use the material after the expiration date.
- Record the LOT number with every process that requires identification of the material.

**The product properties are based on compliance with the instructions for use.**

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**Merz Dental GmbH**