artConcept®

Gebrauchsinformation

- en Instructions for use
- fr Mode d'emploi
- it Istruzioni per l'uso
- es Instrucciones de uso
- pt Instruções de utilização
- pl Instrukcja użycia
- nl Gebruiksaanwijzing
- sv Information för användning
- da Brugsanvisningen
- fi Käyttöohje
- lt Informacinis lapelis
- lv Lietošanas informācija
- et Kasutusteave
- hu Használati tájékoztató
- cs Informace k použití
- sl Navodila za uporabo
- sk Návod na použitie
- el Οδηγίες χρήσης
- tr Kullanım bilgidir

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Fixture and backing of artVeneer® and completion of acrylic teeth, as well as aesthetic correction of dentures.

(en) Instructions for use, please read carefully!

Indication

- Complete loss of individual teeth
- Loss of hard tooth substance

Intended use

Replacement of lost hard tooth substance by:

- Tooth-shaded fixture and for backing artVeneer®
- Veneering telescopic, attachment, bar and implant constructions with artVeneer®
- Tooth-shaded backing of trimmed acrylic teeth
- Fixing acrylic teeth to the model cast with a back protection plate or pontic
- Individualising acrylic teeth, artVeneer® and CAD/CAM restorations produced from M-PM® Disc and the artBloc® Temp
- Repairs

Contraindication

- Hypersensitivity to acrylate compounds

artConcept® consists of

• artPrime PLUS

Metal-adhesive solution that improves adhesion between dental alloys and veneer resins.

• artPreOpaque

Light-curing priming opaquer for use between artPrime PLUS and artOpaque.

• artOpaque

Opaquer, shade-matched to artVeneer® and artDentine to cover metal frameworks, light-curing.

• artDentine/artEnamel Polymer

Tooth-shaded PMMA acrylic, shade-matched to artVeneer® and Merz Dental acrylic teeth.

• artDentine Liquid

Monomer component for artDentine/artEnamel.

• artConnect

Acrylic bonder (artVeneer® - artDentine) especially for highly cross-linked acrylic veneers and PMMA-based acrylic teeth for bonding to artDentine. Also recommended for preparation of mould correction or supplementation of artVeneer®.

• artVeneer®

Multi-layered, tooth-shaded PMMA anterior and posterior veneers, 1:1 combinable with the artegal® tooth lines.

Composition

artPrime PLUS

- Acetone, phosphonic acid monomer, thioctic acid

artDentine Polymer

- PMMA (polymer and copolymer), Barbituric acid catalyst system, organic and inorganic colourants

artDentine Liquid

- MMA (methylmethacrylate), dimethacrylate

artOpaque

- Methacrylic acid ester liquid containing an activator, D,L camphor quinone, Colorants

artPreOpaque

- Urethane acrylates, inorganic fillers, photoinitiators

artConnect

- MMA (methylmethacrylate), THFMA (tétrahydrofurfurylméthacrylate)

Processing

Framework preparation

- Blast the metal frame with 50 - 110 µm aluminium oxide, then rinse with water and dry thoroughly.
- Place a sufficient quantity of **artPrime PLUS** in a mixing bowl; using a disposable brush, apply a layer over the entire adhesive surface and leave for 10 seconds to take effect.

PreOpacuer application

- **artPreOpaque** is applied within 1-3 minutes using a brush with short, stiff bristles (no wash!), especially in the area of the retentions.

Note: No intermediate polymerisation devices may be used for polymerisation.

Opacuer application

- Apply **artOpaque** twice thinly using a brush with short, stiff bristles until the framework is completely covered with colour.

Note: By applying the opaquer thinly, the complete polymerisation of the individual layers is ensured.

Fixture of the veneers

- Remove remnants of wax and fat, grinding dust and other contamination from the **artVeneer®** veneers.
- After evaporation, ensure complete drying, also in the undercuts. There must be no contamination from blowing off with compressed air (e.g. from oil residues in the compressed air system).
- When roughening **artVeneer®** or acrylic teeth, affix small groove retentions as required for larger components.
- Wetten the surface to be joined with **artConnect** and apply sufficient that the connecting surface remains moist for 5 min to work in. If the connecting surface with **artConnect** dries out during this time, moisten again with **artConnect**. The exposure time must not be less than 5 min; the fluid **artDentine** resin has to be applied to the conditioned surface after 10 min at the latest.
- **artEnamel** can be used for individualisation in the incisal area.

Dosage

artDentine/artEnamel

- Mixing ratio artDentine/artEnamel Polymer: Liquid, 1 g : 0.7 ml
- Mixing time 30 sec*
- Dwelling time 15 to max. 45 sec* after mixing
- Fluid phase 240 sec*
- Polymerisation 15 min in the pressure pot at 45 °C and 2-2.5 bar pressure

* Based on a room temperature of 21 °C

Use class IV plaster for making casts where possible. Always sufficiently wetten the models and use standard alginate isolation to separate plaster from acrylic.

Trimming

- Cross-cut carbide burs suitable for resin or suitable cutting discs should be used to rework.
- To prevent inaccurate fitting after polymerisation, avoid excessive heat development during trimming and polishing.

Curing unit artOpaque

DentaColor-XS	Heraeus Kulzer	90 sec
HeraFlash	Heraeus Kulzer	90 sec
Labolight LV-II / III	GC	1 min
LicuLite	Dentsply	1 x HD
Lumamat 100	ivoclar vivadent	Program 2
MPA 200	Dentsply	1 x HD
Solidilite EX	Shofu	1 min
Spektra 2000	Schütz Dental	3 min
Spektra 2000 fast	Schütz Dental	3 min
Spektra LED	Schütz Dental	30 sec
Spektramat	ivoclar vivadent	2 min
Targis Power	ivoclar vivadent	2 min
Triad II	Dentsply	5 min
Uni-XS	Heraeus Kulzer	90 sec

Note

To reduce the residual monomer content, the finished (sanded) but not polished polymerised workpiece should be stored in water at room temperature for about 48 hours before delivery.

Restoration

Products of artDentine can also be restored and supplemented at any time with all standard MMA-based cold-curing polymers. Prior conditioning of the surface is necessary, preferably with artConnect or a similar acrylate-based bonding agent in accordance with the specifications of the manufacturer.

Storage

- Dry and cool area (5 °C to 25 °C)
- Away from direct sunlight

Special notes

- Close all containers tightly after use.
 - Do not use the material after the specified expiration date.
 - Record the LOT number with every process that requires tracking and attribution of the material.
 - If you are using this medical device for a special design, please provide the dentist treating with all the above information.
- Observe the product-specific material safety data sheet with regard to processing.

Risk phrases

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H319 Causes serious eye irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P321 Specific treatment (see on this label).

P362 Take off contaminated clothing.

P501 Dispose of contents/container to.

P210 Keep away from heat. No Smoking.

P233 Keep container tightly closed.

P302+P350 IF ON SKIN: Gently wash with plenty of soap and water.

P271 Use only outdoors or in a well-ventilated area.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.

P403+P233 Store in a well-ventilated place. Keep container tightly closed.

P260 Do not breathe dust/fume/gas/mist/vapours/spray.

P240 Ground and bond container and receiving equipment.

P243 Take action to prevent static discharges.

Disposal

The polymerised components of artConcept® are water-insoluble, inactive, do not pose any threat to groundwater and may therefore be disposed of with normal household waste.

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de.

The product properties are based on compliance with and observation of these instructions for use.

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Merz Dental GmbH