

integral®

& artegral®



ARTIKULATIONS- UND KOMBINATIONSTABELLE

Frontzähne				Seitenzähne	
	Angle I	Angle II	Angle III		
	integral®			integral®	artegral®
OK		UK		OK / UK	OK / UK
CL	UCL	UCM	UCL	(M) L	(M) L (XL)
CM	UCM	UCS	UCL	M	M
CS	UCS	UCS	UCM	S	(XS) S
BL	UBL	UBM	UBL	(M) L	(M) L (XL)
BM	UBM	UBS	UBL	M	M
BS	UBS	UBXS	UBM	S	(XS) S
BXS	UBXS	UBXS	UBS	XS	XS
IL	UIL	UIM	UIL	(M) L	(M) L (XL)
IM	UIM	UIS	UIL	M	M
IS	UIS	UIS	UIM	S	(XS) S