



Frontzahn
Anterior
Anteriore
Antérieure



Formenkarte



Mould Chart

Tabella forme

Ficha de formas

Carte de formes

DeltaForm®




OMP-N® Organic Modified Polymer-Network

Für mehr Qualität For more quality

Per una maggiore qualità


Para una mayor calidad Pour plus de qualité






*OMP-N® specially processed interpenetrated network

Graph based on Ghazal et al.: Two-body wear of resin and ceramic denture teeth in comparison to human enamel. Dental Materials 4 (2007) 24, p. 502-507

Conventional Acrylic Material
Single Cross-Linked Polymer







OMP-N® Organic Modified
Polymer-Network





square






Drei-Schicht-Aufbau
 Three-layer set-up
 Struttura a tre strati
 Estructura de tres capas
 Structure à très couches



16 V-Farben/shades/colori/colores/couleurs A1-D4
 4 Bleach-Farben/shades/colori/colores/couleurs BL1-BL4

Artikulationstabelle DeltaForm®

Cross reference chart · Tabella di articolazione
 Tabla de articulación · Table d'articulation

 DeltaForm®	 DeltaForm®	 DeltaForm® HD
● - Form O1	i6/i7/i8	XL / L
▲ - Form T1 T2	i3/i2 i2/i3/i4/i5	M / S / XS M / S
■ - Form R1 R2 R3 R4 R5 R6	i1/i2/i4 i2/i4/i5 i5/i6/i7 i5/i6/i7 i7/i8 i8	XS / S / M S / M M / L M / L M / L XL
■ - Form S1 S2 S3 S4 S5 S6	i1/i2/i3/i4 i3/i5 i5/i6 i6/i7 i4/i5/i6 i6/i7	XS / S XS / S S / M M / L M / L M / L



Formenkart Seitzenzahn ▶
 Mould Chart posterior
 Tabella forme posteriore
 Ficha de formas posterior
 Carte de formes postérieure
 REF 1090170

Sicherheitsdatenblatt auf Anfrage erhältlich.
 Safety data sheet available on request.
 Fiche de données de sécurité disponible sur demande.
 Scheda dati di sicurezza disponibile su richiesta.
 Puede solicitarse la ficha de datos de seguridad.



Merz Dental GmbH
 Kieferweg 1, 24321 Lütjenburg, Germany
 Tel + 49 (0) 4381/403-0
 Fax + 49 (0) 4381/403-403
 www.merz-dental.de
 EN ISO 13485

