



Seitenzahn
Posterior
Posteriore
Postérieure



Formenkarte


Mould Chart

Tabella forme

Ficha de formas

Carte de formes

DeltaForm®



Lingualisiert

Lingualized

Lingualizzata

Lingualisada

Lingualisé

Hohe Kaustabilität

High masticatory stability

Stabilità masticatoria elevata

Alta estabilidad masticatrice

Haute stabilité masticatrice


Geringe Horizontalschübe

Less horizontal forces

Spinte orizzontali minori

Fuerzas horizontales reducidas

Forces horizontales réduites




Mörser-Pistill-Prinzip - Selbstregulierende Okklusion - Einfache Logik des Einschleifens

The mortar pestil principle (cupola-dip) - Self-regulating occlusion - Simple logic of grinding

Principio „mortajo-pestello“ (superfici concave-convexe) - Occlusione autoregolante - La logica semplice della fresatura

Principio del mortero y la mano (cúpula y concavidad) - Oclusión auto regulada - Lógica sencilla en el ajuste oclusal

Principe du mortier-pilon (cuspidé-fosse) - Occlusion à autorégulation - Logique simple du meulage




Moderate, geringe Höckerneigung von 20°-15°, distal abnehmend

Moderate, little cusp angle of 20°-15°, towards distal declining

Inclinazione moderata delle cuspidi, di 20°-15°, decrescendo distalmente

Moderada, ligera inclinación cuspidéa de 20°-15°, decreciente hacia distal

Inclinaison modérée faible des cuspides de 20°-15°, décroissance distal




Patentierte Abrasionsfacetten, altersgerechte Kontur, geringes Einschleifen

Patent abrasion facets, age-based contour, less grinding

Facette abrasive brevettate, contorni in funzione all'età, molaggio ridotto

Facetas de abrasión patentadas, perfil adecuado a la edad del paciente, requieren poco tallado

Facettes d'abrasion brevetées, contour adapté à l'âge du patient, faible meulage



Patentierter Preventive Shift Groove (PSG[®]) - An den Kieferkamm angepasste basale Kontur

Patent Preventive Shift Groove (PSG[®]) - To the ridge fitted basal contour

Preventive Shift Groove (PSG[®]) brevettato - Contorno basale adattato alla cresta mascellare

Sistema patentado Preventive Shift Groove (PSG[®]) - Perfil basal adaptado a la cresta alveolar

Preventive Shift Groove (PSG[®]) brevetée - Contour basal adapté à la crête








- produzione
A punctero

- fabricación
Mise au point

- fabrication



DeltaForm® / DeltaForm® HD






Drei-Schicht-Aufbau


Three-layer set-up
Structura a tre strati
Estructura de tres capas
Structure à très couches

16 V-Farben/shades/colori/colores/couleurs A1 - D4
4 Bleach-Farben/shades/colori/colores/couleurs BL1 - BL4

DeltaForm®

				
inTEGRAL®	arTEGRAL®/arTEGRAL® Life	inTEGRAL®	arTEGRAL®/arTEGRAL® Life	DeltaForm® / HD
CL CM CS	CL CM CS	UCL UCM UCS	UCL UCM UCS	XL/L M S/XS
BL BM BS	BXL BL BM BS BOL BOM BOS BIL BIM BIS	UBLL UBL UBM UBS	UBXL UBL UBM UBS UBIL UBIM UBIS UBIL UBIM UBIS	XL XL/L M S/XS XL/L M S/XS XL/L M S
IL IM IS	IXL IL IM IS ITL ITM ITS	UIL UIM UIS	UBXL UIL UIM UIS UIL UIM UIS	XL L M S/XS L M S/XS
Polystar® EDITION <i>Selection</i>		Polystar® EDITION <i>Selection</i>		DeltaForm® / HD
CL2 CM2 CM4 CM6 CS2		UL2/UL4 UM2/UM4 UM2/UM4 UM6/US2 US2		XL/L L/M M/S M/S S/XS
BL1 BL2 BM2 BM4 BM6 BS2		UL2 UL2/UL4/UM4 UM2/UM4 UM2/UM4 UM6/US2 US2		XL/L XL/L L/M M/S S/XS XS
BIL2 BIM2 BIM4 BIM6 BIS2		UL2/UL4 UM2/UM4 UM2/UM4 UM6/US2 US2		M/S M/S M/S S/XS XS
IL2 IM2 IM4 IM6 IS2		UL2/UL2/UL4 UM2/UM4 UM2/UM4 UM6/US2 US2		L/M L/M M/S M/S S/XS

Formenkarte Frontzahn ▶
Mould Chart anterior
Tabella forme anteriore
Ficha de formas anterior
Carte de formes antérieure
REF 1090171



Merz Dental GmbH
Kieferweg 1 24321 Lütjenburg, Germany
Tel + 49 (0) 4381/403-0
Fax + 49 (0) 4381/403-403
www.merz-dental.de
EN ISO 13485

CE 0482

